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**GB**  
Ureteral balloon dilation catheter  
Instructions for use

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**IFU Nr. 29-A**

Symbols see part B

**Instructions for use**  
Ureteral balloon dilation catheter

IFU Nr 29-A Dilatationskatheter EN\_02/17.02.2016

Please note: This medical product may be obtained and used only by trained medical personnel.

**1 Description**  
Ureteral balloon dilation catheter

**2 Contents and packaging**  
- 2-way balloon dilation catheter

Figure 1: Hydrophilic part of the balloon dilation catheter

- Instructions for use

**3 Intended use**  
Dilatation of the ureter.

**4 Indications**  
Dilatation of ureteral strictures or for ureteral dilatation before ureteroscopy or treatment of stones.

**5 Contraindications**  
Contraindications result from restrictions for the ureteroscopy, for example: Serious infections, urosepsis, prolonged obstructions, chronic urinary obstruction or extrinsic obstructions.

**6 Side effects**  
- Hematuria  
- Foreign body irritation in the bladder and renal regions  
- Erosions and ulcerations of the ureter  
- Dysuria  
- Subsequent urge to urinate  
- Feeling of pressure in the flank  
- Urinary retention (caused by ostial swelling)  
- Urinary tract infection  
- Sepsis

**7 Instructions**  
**a. Preparing the balloon dilation catheter**  
Before use, the catheter must be inspected to see if it is intact (packaging, product).  
Make sure that all air is removed from the balloon before use:  
1. Carefully remove the balloon's protective cover.  
2. Fill the balloon with approx. 4 ml saline solution.  
3. Empty the balloon completely using a syringe. Repeat the filling and emptying procedure until the balloon is completely void of air.  
4. Activate the hydrophilic coating with sufficient moistening from saline solution or sterile water.

**b. Insertion/removal**  
1. Lay a guide wire through to the renal pelvis using radiological control.  
2. Carefully advance the balloon dilation catheter over the guide wire (note radiopaque markers).  
3. After placing the balloon dilation catheter, replace the syringe with a high-pressure syringe filled with saline solution and contrast agent.  
4. Completely empty the balloon before removal.

**8 Possible complications and/or risks**  
- Exceeding the maximum allowable balloon inflation pressure may lead to tissue trauma, damage to the ureter wall, rupture or mucosa lesions.  
- Exceeding the Rated Burst Pressure (RBP) may damage the balloon and/or catheter (see table).  
- The balloon size is selected by the surgeon depending on the ureteral lumen.  
- The dilatation time should not exceed 5-10 minutes in order to prevent ischemia of the surrounding tissue.

**9 Instructions for use / precautions**  
- This product is only sterile if the packaging is undamaged and unopened.  
- For single use only!  
- Do not use products after the date of expiry has passed.  
- Do not wipe the product with alcohol/disinfectant agents. This may damage the surface.  
- All components must be carefully checked for compatibility and integrity before use.

- Do not damage the products with sharp surgical instruments – this increases the risk of catheter breakage or damage to the balloon.
- 5 FR introducer sheaths are endoscope compatible with 4.0 and 6.0 mm balloon diameters, and 6 FR sheaths are compatible with an 8 mm balloon diameter.
- The correct position of the catheter and the complete expansion of the balloon must be radiologically verified.
- The hydrophilic properties of the balloon dilation catheter are only activated if the catheter is sufficiently moistened with physiological saline solution before use. Failure to do this will make it extremely difficult to push the guide wire forwards.
- The balloon and lumens for inflation and the guide wire must be void of air.
- Only fill the balloon with a mixture of saline solution and contrast agent (1:1). Never fill the balloon with air or other gases.
- For guide wires in situ, check that the guide wire diameter is compatible with the specified inner lumen of the catheter (see label)
- Do not move the guide wire while the balloon is being filled.
- Lack of proper care when advancing the catheter may damage the ureter.
- In exceptional cases, it may be impossible to unblock the balloon. In this case, the balloon can be emptied by cutting off the valve.
- The balloon may only be inserted and removed when it is empty.

**10 Interaction with other drugs**  
The catheter could potentially discolor due to the concomitant administration of medication containing dye-based urine presentation media, or even as a result of contact with urine.

**11 Transport and storage conditions**  
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport.  
Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

**12 Disposal**  
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

**Overview of the pressure tolerances:**

Inflationsdruck		Ballondurchmesser (mm)		
atm / bar	kPA	4.0	6.0	8.0
4	400	3.78	5.68	7.45
6	600	4.00	6.00	8.00
8	800	4.15	6.16	8.31
10	1000	4.28	6.30	8.53
11	1100	4.33	6.37	8.64
12	1200	4.40	6.44	8.77
13	1300	4.46	6.52	8.89
14	1400	4.53	6.59	9.00
15	1500	4.60	6.67	9.12
16	1600	4.67	6.79	9.25
17	1700	4.75	6.84	

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