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GB  
Ureteral balloon occlusion catheter  
Instructions for use

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## IFU Nr. 34-A

Symbols see part B

### Instructions for use Ureteral balloon occlusion catheter

Please note: These medical products may be obtained and used only by trained medical personnel.

#### 1. Description

Ureteral balloon occlusion catheter with detachable connector.

#### 2. Contents and packaging

- Double-lumen ureteral balloon occlusion catheter with hydrophilic tip, integrated mid-balloon x-ray marker and detachable connector.

- Protective tube

- Disposable packaging with instructions for use

#### 3. Intended use

The balloon occlusion catheter is used for temporary closure of the proximal ureter or pyeloureteral junction and for injecting solutions, e.g. contrast medium.

#### 4. Indication

##### Therapeutic indication:

- Percutaneous nephrolithotripsy

##### Diagnostic indication:

- Application of contrast medium

#### 5. Contraindications

Contraindications result from restrictions for the ureterorenoscopy, for example: Serious infections, urosepsis, prolonged obstructions, chronic urinary obstruction or extrinsic obstructions.

#### 6. Possible side effects

- Hematuria  
- Foreign body irritation in bladder and renal regions  
- Erosions and ulcerations to the ureter  
- Dysuria  
- Subsequent urge to urinate  
- Feeling of pressure in the flank  
- Urinary retention (caused by ostial swelling)  
- Urinary tract infection  
- Sepsis

#### 7. Instructions

The balloon occlusion catheter should be introduced using endourological standard method (with direct radiological fluoroscopy).

#### Assembly and use:

1. Remove the protective tube and activate the hydrophilic coating at the catheter tip with sufficient moistening from saline solution or sterile water.

**Remove the stylet prior to use if the guide wire is still in situ (loosen the yellow cap and remove stylet). If too much resistance is encountered when advancing via the guide wire, the connector can be loosened (see Figure 1). It must be retightened once the catheter has been positioned! (see Figure 3)**

2. Introduce the balloon occlusion catheter via the working channel of a positioned cystoscope (the correct positioning of the balloon can be checked using x-ray to detect the marker (at balloon midpoint)).

3. After the desired position has been achieved, loosen the connector (see Figure 1). Then remove the connector, including the stylet


Figure 1: Loosening the connector

4. Removing the cystoscope  
To inflate the balloon, place the connector back on the end of the catheter (until it comes to a stop, see Figure 2).


Figure 2: Pushing the connector until it comes to a stop

6. Fixing the connector in place (see Figure 3)


Figure 3: Fixing the connector in place

- Before filling the balloon, set the empty syringe at the filling vent and create a vacuum in the balloon. Then take the syringe away. This will prevent air from entering the balloon as it is being filled.
- Set the syringe at the filling vent and block the balloon with 1 ml of NaCl 0.9% with x-ray monitoring.
- Inject contrast medium via the Luer lock.

- Removal of the occlusion catheter:**
  - Passive emptying by removing the connector.
  - Active emptying by attaching an empty syringe to the filling vent and carefully aspirating the blocking fluid.
  - The balloon occlusion catheter can be removed after deblocking has been completed.

- Possible complications and/or risks**
  - Blocking for too long can lead to ureter tissue damage.
  - The catheter cannot be pushed through the ureter and into the kidney due to a blockage.
  - Rough-edged stone deposits may cause the balloon to burst.
  - Injury to the urethra or ureter resulting in a subsequent stricture
  - Injury to the bladder

- Warnings / Precautions**
  - This product is only sterile if the packaging is undamaged and unopened.
  - For single use only!
  - Do not use catheter after the date of expiry has passed.
  - Do not wipe the product and set components with alcohol/disinfectant agents. This may damage the surface.
  - All components must be carefully checked for compatibility and integrity before use.
  - The connector should not be screwed too tight in order to prevent damage to the tube.
  - Lack of proper care when advancing the balloon occlusion catheter may damage the ureter.
  - Do not damage the catheter or stylet with sharp surgical instruments - this causes increased risk of breakage.
  - Do not excessively stretch the catheter before use. Stretching it too far increases the risk of breakage.
  - Catheter and stylet must not be sharply curved (risk of snapping).
  - Check that the guide wire diameter is compatible with the stated inside lumen of the catheter for guide wires in situ.
  - After repositioning the patient check the balloon's position and blocking.
  - Do not block the balloon with contrast medium.
  - Do not over-block the balloon! Observe the maximum volume listed on the label.
  - The balloon occlusion catheter is intended for a temporary service life of ≤ 60 minutes.

- 10. Interaction with other drugs**

The catheter could potentially discolor due to the concomitant administration of medication containing dye-based urine presentation media, or even as a result of contact with urine.

- 11. Transport and storage conditions**

The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport.

Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

- 12. Disposal**

After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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