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"Tiger" guide wire with PTFE coating "Semislip" with hydrophilic tip, "Visiomulti" guide wire with 2 flexible ends
Instructions for use

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IFU Nr. 19-A

Symbols see part B

Instructions for use "Tigerdraht"™ guide wire with PTFE coating, "Semislip"™ with hydrophilic tip, "Visiomulti"™ guide wire with 2 flexible ends

IFU Nr 19-A UTD_UAS Guide Wire_EN_02/06.05.2015

Please note: This medical product may be obtained and used only by trained medical personnel.

1. Description
Guide wires are used for all instruments and catheters that are inserted in a minimally-invasive procedure.

2. Contents and packaging
- Guide wire made of nitinol with yellow-black coiled PTFE coating, 0.035 inches (Tigerdraht).
- Guide wire made of nitinol with PTFE coating, flexible tip 5 cm made of Pebax, hydrophilically coated (semislip).
- Guide wire made of nitinol with white-black coiled PTFE coating 0.035 inches, super-elastic with 2 flexible ends (visiomulti).

3. Intended use
A guide wire facilitates the placement of a catheter or instrument during diagnostic or therapeutic applications in the urinary tract.

4. Indications
Endourological applications: (Ureterscopy, ureteral stenting, insertion of a ureteral or nephrostomy catheter, etc.)

5. Contraindications
Contraindications arise from restrictions on the specific applications

6. Side effects
- Hematuria
- Infection (ascending)
- Hemorrhages
- Bacteremia
- Dysuria
- Flank pain
Additional side effects arise from individual specific applications.

7. Instructions
Remove the guide wire with the dispenser tube from the packaging, then moisten the hydrophilic guide wire and remove it from the dispenser tube.

Internal ureteral stenting:
Ureteral stenting one-sided closed
- Push the ureteral stent with the open end onto the rigid end of the guide wire.
- Ureteral stent is attached to the guide wire by a clamp.
- Then insert via the working channel of the cystoscope and advance up to the renal pelvis.
- After placing the stent, the wire is removed.

Ureteral stent open on both sides
- Guide the guide wire, flexible tip first, up to the renal pelvis.
- Place the ureteral stent using the guide wire.
- Removal of the guide wire.

Ureteral stent intraoperatively
- Open the ureter at the level of the renal pelvis and push the catheter into the bladder using the guide wire.
- Remove the guide wire and then insert it through the drainage hole in the middle of the stent.
- Thereby, place the stent in the kidney and remove the guide wire.

Percutaneous nephrostomy:
- Initial puncture and insertion of the flexible end of the guide wire.
- Then, the guide wire is placed in the renal pelvis.
- Attachment of the wire and careful removal of the puncture needle.
- The dilator/catheter is then inserted using the placed guide wire.

8. Possible complications and/or risks
- Damage to the urethra or ureter resulting in a subsequent stricture.
- Damage to the bladder
- Damage to the kidney due to uncontrolled advancement into the renal pelvis.
- The guide wire cannot be pushed through the ureter into the kidney due to a blockage.
- Ureter spasms
- Urinary tract infection
- Sepsis

9. Warnings / Precautions
- This product is only sterile if the packaging is undamaged and unopened. For single use only!
- Observe shelf life (see label).
- Do not clean product and set components with alcohol or treat with disinfectant. This may damage the surface.
- Do not re-sterilize! This may cause material fatigue and inadequate sterility!
- Compatibility with products from external manufacturers (do not take set components into account).
- All components must be carefully checked before use. Products that may be damaged must not be used.
- This guide wire should only be used by medical personnel who have been trained to use guide wires.
- The guide wire must not be withdrawn or manipulated over sharp edges, such as those found in cannulas, as this can lead to the destruction and/or dissolution of the external Teflon cladding.
- Do not use a wire diameter larger than the specified internal lumen of the stent, catheter, puncture cannula or dilator.
- The guide wire may not be pushed rigid end first, to prevent traumatization of the renal pelvis or bladder wall.
- In the case of a possible deadlock of the wire with the stent, the wire must be corrected immediately under all circumstances.
- A sufficient guide wire end must remain hanging out to ensure good access to the guide wire at all times.
- Only apply the guide wire under visual guidance, to monitor the position of the guide wire.
- The guide wire can be damaged by uncontrolled pushing or pulling over rigid, sharp edges, thereby stripping away the plastic coating.
- The hydrophilic wires must be moistened prior to use. A spraying adapter is installed on the dispenser tube for this purpose.
- Hydrophilic wires cannot be allowed to dry out during use, as this may lead to sticking and thus damage to the instrument - please take into account if application duration is prolonged.

10. Interaction with other drugs
The guide wire could potentially be discolored when applying urinary examination methods that use dye.

11. Transport and storage conditions
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

12. Disposal
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations.

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