DE

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Guide wires

IFU Nr. 26-A Führungsdrähte EN/02/05.05.2015

Please note: This medical product may be obtained and used only by trained medical personnel.

1. Description
Guide wires are used for all instruments and catheters that are inserted in a minimally-invasive procedure.

2. Contents and packaging (depending on order number)
Guide wire PTFE coated (Lunderquist guide wire)
Guide wire PTFE coated (mobile core)
Guide wire PTFE coated (standard)

3. Intended use
A guide wire facilitates the placement of a catheter or instrument during diagnostic and therapeutic applications in the urinary tract.

4. Indications
Endourological applications: Ureterscopy, ureteral stenting, insertion of a ureteral or nephrostomy catheter, or changing of a suprapubic catheter, etc.

5. Contraindications
Contraindications arise from restrictions on the specific applications

6. Side effects
- Hematuria
- Infection (ascending)
- Hemorrhages
- Bacteremia
- Dysuria
- Flank pain
Additional side effects arise from individual specific applications.

7. Instructions
Internal ureteral stenting:
Ureteral stenting one-sided closed
- Push the ureteral stent with the open end onto the rigid end of the guide wire.
- Ureteral stent is attached to the guide wire by a clamp.
- Then insert via the working channel of the cystoscope and push forward up to the renal pelvis.
- After placing the stent, the wire is removed.
Ureteral stent open on both sides
- Guide the guide wire, flexible tip first, up to the renal pelvis.
- Place the ureteral stent using the guide wire.
- Removal of the guide wire.
Ureteral stent intraoperatively
- Open the ureter at the level of the renal pelvis and push the catheter using the guide wire in the bladder.
- Remove the guide wire and then insert it through the drainage hole in the middle of the stent.
- Thereby, place the stent in the kidney and remove the guide wire.
Percutaneous nephrostomy:
- Initial puncture and insertion of the flexible end of the guide wire.
- Then, the guide wire is placed in the renal pelvis.
- Attachment of the wire and careful removal of the puncture needle.
- The dilator/catheter is then inserted using the guide wire.
Suprapubic catheter change:
- Insert the exchange wire, flexible tip first, into the bladder. If the wire can't be easily inserted-->hold the catheter as straight as possible (without bends).
- Remove the catheter from the bladder using the placed exchange wire.
- When using dilation exchange sets, the wire can be used for dilation.
- Insert the catheter into the bladder using the placed exchange wire.
- After placing the catheter, the wire is withdrawn and removed.

8. Possible complications and/or risks
- Damage to the urethra or ureter resulting in a subsequent stricture.
- Damage to the bladder
- Damage to the kidney due to uncontrolled advancement into the renal pelvis.
- The guide wire cannot be pushed through the ureter into the kidney due to a blockage.
- Ureter spasms
- Urinary tract infection
- Sepsis

9. Warnings / Precautions
- This product is only sterile if the packaging is undamaged and unopened. For single use only!
- Observe shelf life (see label).
- Do not clean product and set components with alcohol or treat with disinfectant. This may damage the surface.
- Do not re-sterilize! This may cause material fatigue and inadequate sterility!
- Compatibility with products from external manufacturers (do not take set components into account).
- All components must be carefully checked before use. Products that may be damaged must not be used.
- This guide wire should only be used by medical personnel who have been trained to use guide wires.
- The guide wire must not be withdrawn or manipulated over sharp edges, such as those found in cannulas, as this can lead to the destruction and/or dissolution of the external Teflon cladding.
- Do not use a wire diameter larger than the specified internal lumen of the stent, catheter, puncture cannula, or dilator.
- The guide wire may not be pushed rigid end first, to prevent traumatization of the renal pelvis or bladder wall.
- In the case of a possible deadlock of the wire with the stent, the wire must be corrected immediately under all circumstances.
- A sufficient guide wire end must remain hanging out to ensure good access to the guide wire at all times.
- Only apply the guide wire under visual guidance, to monitor the position of the guide wire.
- If the wire tip is in a curved vessel form, the movable core must not be pushed forward.
- The core must never be rotated or pushed forward with force, as too much pressure could puncture the wire coil.

10. Interaction with other drugs
The guide wire could potentially be discolored when applying urinary examination methods that use dyes.

11. Transport and storage conditions
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

12. Disposal
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable regulations and guidelines.

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