DE
UROLIX Hydrophilic Guide Wires
GebrauchsweisungEN
UROLIX Hydrophilic Guide Wires
Instructions for useFR
Fils-guides hydrophiles UROLIX
Mode d'emploiES
Alambres guia hidrofilicos UROLIX
InstruccionesIT
Filii guida idrofili UROLIX
IstruzioniPT
UROLIX Fio-guias hidrofilicos
Instruções de utilizaçãoDA
UROLIX Hydrofile guidewires
Betjeningsvejledning

IFU Nr. 36-A

EN
Instructions for use
UROLIX Hydrophilic Guide WiresThis IFU can also be downloaded from our website: urotech.com
IFU Nr. 36-A Urolinx_EN_01/14.06.2019

Please note: This medical product may be obtained and used only by trained medical personnel. Non-intended use of the product can cause serious and irreversible damage.

1. Description
Guide wires are used for all instruments and catheters that are inserted in an endo-urological procedure.

2. Contents and packaging
Guide wire made of nitinol with hydrophilic coating

3. Intended use
A guide wire facilitates the placement of a catheter or instrument during diagnostic or therapeutic applications in the urinary tract.

4. Indications
Endourological applications: (Ureteroscopy, ureteral stenting, insertion of an Ureteral Access Sheath or Ureterorenoscope, etc.)

5. Contraindications
Contraindications arise from restrictions for an endo-urological procedure:

- Urinary Tract Infection
- Urosepsis
- Uncorrect bleeding diastasis
- Prolonged obstructions, chronic urinary obstruction, extrinsic obstructions

6. Patient population

The patient population or target patient group is derived from the indication determined by the responsible doctor who therapeutically treats the patient as part of an endo-urological procedure (the leading procedure per se) according to the intended use of the medical device. There are no known restrictions to the patient population or the target patient group.

Use of the product on minors:

The user can use the product on minors in the case that the physiological and anatomical conditions of the patient permit the use of the product.

Use of the product on women who are pregnant or breastfeeding:
The indication for the use of the product on women who are pregnant or breastfeeding must be narrowed down by the user based on the respective individual physiological and anatomical conditions of the patient.

7. Possible side effects

- Hematuria
- Infection (ascendant)
- Hemorrhages
- Bacteremia
- Dysuria
- Flank pain

Additional side effects arise from individual specific applications.

8. Instructions

Remove the guide wire with the dispenser tube from the packaging, then moisten the hydrophilic guide wire and remove it from the dispenser tube.

Internal ureteral stenting:

Ureteral stenting one-sided closed

- Push the ureteral stent with the open end onto the rigid end of the guide wire.
- Ureretal stent is attached to the guide wire by a clamp.
- Then insert via the working channel of the cystoscope and advance up to the renal pelvis.
- After placing the stent, the wire is removed.

Ureteral stent open on both sides

- Guide the guide wire, flexible tip first, up to the renal pelvis.
- Place the ureteral stent using the guide wire.
- Removal of the guide wire.

Insertion of Access Sheath_Ureterorenoscope etc.

- Get access to the Orifice in the bladder by using a cystoscope.
- Insert the Guide Wire through the cystoscope into the Ureter and the kidney, using X-ray for verification.
- Remove the cystoscope
- Insert the Access Sheath, Ureterorenoscope etc. over the Guide Wire into the upper urinary tract.

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Urotech GmbH | Medi-Globe-Straße 1-5