DE

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
GB

**Instructions for use**  
**Ureteral access sheath**


IFU Nr 35-A Zugangsschleuse EN\_02/23.02.2016

*Please note: This medical product may be obtained and used only by trained medical personnel.*

- Description**  
Ureteral access sheath
- Contents and packaging**  
- Ureteral access sheath with dilator  
- Instructions for use
- Intended use**  
The ureteral access sheath is indicated for use in endoscopic urological procedures and facilitates the passage of endoscopes and other instruments through the urinary tract.
- Indication**  
Diagnostic/Therapeutic flexible ureterorenoscopy (URS)
- Contraindications**  
Contraindications result from the restrictions for the ureterorenoscopy.
- Possible side effects**  
Using the access sheath may cause mucous membrane irritation or inflammations in the urinary tract in some circumstances.
- Instructions**  
a) Guide the dilator into the sheath and slide the dilator handle onto the sheath funnel as illustrated. The dilator is securely fixed for use as soon as the red indicator appears in the viewing window. Check the tip transition to the sheath tube. This should be atraumatic (see illustration).

- Wet the sheath with saline solution or sterile water to activate the hydrophilic coating.
- Guide the sheath unit carefully over the inserted guide wire and into the ureter.
- Now remove the dilator by pushing the two rear fins of the dilator unit together as illustrated and then pulling it out of the sheath.

- Now guide the necessary instruments through the access sheath depending on requirements.
- If a retrograde pyelogram is required, the dilator can be reintroduced and then contrast medium injected through the Luer connection on the dilator handle.

- Possible complications and/or risks**  
- Kidney stones may be too large for removal using an access sheath and may get stuck in the tube as a result.  
- Careless protrusion of the guide wire or the access sheath may cause damage to the urinary tract as well as perforation of the kidney, the renal pelvis, the ureter or the bladder.  
- If the sheath is dislocated, improper repositioning (see precautions) may cause significant damage to the ureter.
- Warnings / Precautions**  
- This product is only sterile if the packaging is undamaged and unopened.  
- For single use only!  
- Do not use the product after the expiry date has passed.  
- Do not wipe the product and set components with alcohol/disinfectant agents. This may damage the surface.  
- The compatibility and integrity of all of the components and parts required for the procedure should be checked carefully prior to use.  
- If the access sheath is dislocated, it must be completely repositioned! To prevent damage, the sheath must be completely removed and combined with the dilator unit again from outside of the body.  
- Do not forcibly remove the dilator from the sheath! If it is not possible to pull the dilator unit out of the sheath, remove the entire product.  
- Do not clamp the access sheath; use a sealing plug if necessary.  
- Expanding the ureter and urethra more than necessary causes bleeding at the site of application and may lead to tissue damage.  
- Use only saline solution or sterile water to activate the hydrophilic coating.  
- Introduce instruments with caution to avoid damaging the sheath.  
- Confirm the correct position of the sheath using X-ray.  
- The ureter sheath is indicated for an application of less than 60 minutes.

**10. Cross reactions**  
It is possible that the sheath may become discolored due to the use of dye-based urine imaging methods, or even as a result of contact with urine.

**11. Transport and storage conditions**  
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

**12. Disposal**  
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.