EN

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EN

Instructions for Use Renodrain and Magnusson nephrostomy direct puncture set

IFU Nr.14-D Direktpunktion EN 08/12.04.2017

Please note: This medical product may be obtained and used only by trained medical personnel.

1. Description

Renodrain and Magnusson locking puncture sets with pigtail catheter for percutaneous nephrostomy.

- 2. Contents and packaging**
- Nephrostomy catheter
 - Schüller guide wire
 - Mandrin or obturator
 - Universal adapter
 - Catheter plug
 - 2 or 3 part puncture needles
 - Double packaging and instructions for use

3. Intended use

Percutaneous urinary diversion from the kidney.

4. Indications

Note: Before determining the indication for percutaneous nephrostomy, it should first be established whether the urinary obstruction can be adequately managed using an internal ureteral stent.

- Diagnostic Indications:**
- Antegrade pyelography and ureteroscopy
 - Collection of material for investigation (histology, cytology, microbiology)
 - Urodynamic measurements in the upper urinary tract
 - Determination of renal function
- Therapeutic Indications:**
- Decompression in cases of obstructive uropathy of the upper urinary tract
 - Urinary diversion in the case of a ureteral fistula
 - Initial step towards antegrade manipulations in the upper urinary tract (ureteral splinting, ureteral embolization, PNL, laser treatment, ureterotomy)

5. Contraindications

- Coagulation disorders
- Severe inflammation
- Dermatological conditions around the puncture site
- Malformed kidneys
- Severe spinal curvature

Contraindications result from restrictions for the corresponding basic applications.

6. Possible side effects

Bleeding, vascular lesions, loss of parenchyma, hematuria, infection (ascending), pain at the puncture site, pain in connection with the lack of urine production, perinephric extravasation and bacteremia.

7. Instructions

When using this product, the general hygiene rules applicable to catheterization must be complied with.

Preparation, storage and patient anesthesia:

- Controlled conditions in the ward, sterile conditions with washing and sterile covering as for open surgery or PNL
- Patient in abdominal position on the table
- Antegrade PNS, generally in local infiltration anesthesia or, in the case of a retrograde PNS or extension of the intervention to a percutaneous operation, regional anesthesia or intubation anesthesia, depending on the patient's resilience and the surgeon's experience.

Technical execution:

Instructions (preparation for one-step technique)

Note: The sharpened tip of the cannula may not be used to stretch the pigtail as this may damage or puncture the inside of the catheter.

1. Stretch the pigtail with the stretching sleeve located on the catheter
2. Replace the inner part of the cannula with the mandrin and insert it into the catheter until it comes into contact with the Luer-lock connector
3. Screw the cannula connection together
4. Remove the stretching sleeve through the tip of the catheter
5. Remove the mandrin once again and replace it with the inner part of the cannula.
6. The catheter is now ready for use

Technical execution of the one-step technique:

- Localization of the kidney / renal pelvis under ultrasound or X-ray guidance.
- Precise initial puncture at the posterior axillary line below the costal arch in the direction of the lower lateral renal circumference.
- After the puncture, check if the renal pelvis has in fact been punctured by removing the inner part of the cannula (indicator: urine runs out).
- If otherwise, insert the needle into the cannula once more and pierce again.
- After the renal pelvis has been successfully punctured, pull the cannula out of the catheter and allow the catheter to drain. The catheter assumes its predetermined pigtail shape when the cannula is removed.
- Attach the catheter to the skin with sutures.

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Instructions (two-step technique)

1. Stretch the pigtail with the stretching sleeve located on the catheter
2. Insert the obturator, remove the mandrin, screw the Luer-lock connector together
3. Remove the stretching sleeve through the tip of the catheter
4. The catheter is now ready for use

Technical execution of the two-step technique:

- Localization of the kidney / renal pelvis under ultrasound or X-ray guidance.
- Precise initial puncture at the posterior axillary line below the costal arch in the direction of the lower lateral renal circumference.
- Insert the guide wire, flexible tip first, up to the renal pelvis.
- Positioning the drainage catheter using an indwelling wire.
- Checking the position of the catheter by imaging with contrast medium.
- Removing the guide wire and mandrin once exact positioning has been achieved.
- The distal catheter end assumes its predetermined pigtail shape. Securing the catheter with skin suture.

Magnusson Locking Pigtail

Note: For puncturing, see the two-step technique

- After placing the catheter, draw the wire back slowly and at the same time pull the movable connector of the Magnusson locking pigtail system in the direction of the fixed connector and twist to lock these together. In this way, the pigtail is fixed securely in the kidney.
- Note: To securely attach the free end of the suture to the catheter, it is recommended to turn the movable connector by 360°.
- To release the pigtail, turn the connection further to the right and push the movable connector towards the body again.

8. Possible complications and/or risks

- Puncturing blood vessels in the wall of the kidney can lead to hematuria.
- Incorrect puncture
- Incrustation, blood clotting, and kinks in the catheter can lead to partial or total obstruction of the suprapubic catheter with the risk of infected hydronephrosis, and in extreme cases urosepsis.
- Dislocation of the catheter.
- Leakage of urine via the nephrostomy channel after removal of the catheter.
- With the one-step technique, correcting the position of the catheter means performing a new puncture.
- Risk of perforation of the kidney and the renal pelvis.

9. Instructions for use / precautions

- This product is only sterile if the packaging is undamaged and unopened.
- For single use only! Do not re-use, recycle, or re-sterilize. Re-using, recycling, or re-sterilizing can impair the structural integrity of the instrument and/or lead to malfunctioning of the instrument, which can in turn lead to illness, injury, or death of the patient.
- All components must be carefully checked for compatibility and integrity before use.
- Do not use products after the date of expiry has passed.
- Do not wipe the product and the set components with alcohol/disinfectant agents. This may damage the surface.
- Do not damage the catheter with sharp surgical instruments - this increases the risk of breakage.
- The guide wire can be damaged by uncontrolled pushing or pulling over rigid, sharp edges, thereby stripping away the plastic coating.
- The sharpened tip of the cannula (one-step technique) may not be used to stretch the pigtail as this may damage or puncture the inside of the catheter.
- The needle can be bent due to excessive force while puncturing. The needle can be stabilized by supporting or fixing it in the puncture area, thereby able to prevent bending or kinking.
- Check the catheter regularly for drainage effectiveness and position.
- Only lubricate the catheter with water-soluble lubricants. Do not use petroleum-based lubricants.
- Catheters must be changed after no more than 30 days.

10. Interaction with other drugs

The catheter could potentially discolor due to the concomitant administration of medication containing dye-based urine presentation media, or even as a result of contact with urine.

11. Transport and storage conditions

The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport.

Products must be stored dry and protected from direct sunlight, in a temperature range of 5 - 30 °C.

12. Disposal

After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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