	DE	GB	FR	IT	ES	PT	GR	AE	CN	JP	CZ
	Hersteller	Manufacturer	Fabricant	Fabbricante	Fabricante	Fabricante	Κατασκευαστής	الجهة الصانعة	製造商	メーカー	Výroce
	Herstellungdatum	Date of manufacture	Date de fabrication	Data di fabbricazione	Fecha de fabricación	Data de fabricaço	Ημερομηνία παραγωγής	تاريخ التصنيع	生产日期	製造日	Datum výroby
	EU-Konformitätszeichen (Nummer gibt die Registrierungsnummer der benannten Stelle an)	EU conformity symbol (number indicates the registration number of the notified body)	Marquage CE de conformité (les chiffres indiquent le numéro d'identification de l'organisme notifié)	Marchio di conformità UE (il numero indica il numero di registrazione dell'organismo notificato)	Marca de conformidad de la UE (el número indica el número de registro de organismo designado)	Símbolo de conformidade com a legislação europeia (indica o número de registo do organismo notificado)	Σύμβολο συμμόρφωσης ΕΕ (ο αριθμός είναι ενδεικτικός για τον αριθμό καταχώρησης του κοινωπονημένου φορέα)	رمز مطابقة الاتحاد الأوروبي (الأرقام إلى رقم تسجيل الهيكل المعنى)	歐盟合格標誌 (數字所指相 关认证机构的 注册号)	CEマーキング (番号は通知機関の登録番号を示す)	EU značka shody (číslo je číslo registrátním číslém jmenovaného nisa)
	Bestell-nummer	Order number	Numéro de commande	Numero d'ordinazione	Nº. de pedido	Numero do pedido	Αριθμός παραγγελίας	رقم طلب الشراء	订单号	注文番号	Objednací číslo
	Chargen-bezeichnung	Batch code	Numéro de lot	Numero del lotto	Código de lote	Numero de lote	Χαρακτηρισμός περιθώ	اسم الطلبيّة	批号	バッチコード	Označení šarže
	Verwendbar bis	Use by	Utilisable jusqu'au	Utilizzabile entro	Fecha de caducidad	Prazo de validade	Ημερομηνία λήξης	صالحه للاستخدام حتى	有效日期	有効日期	Použitelné do
	Sterilisation mit Ethyloxid	Sterilisation with ethylene oxide	Sterilisation par oxyde d'éthylène	Sterilizzazione con ossido di etilene	Esterilización con óxido de etileno	Esterilização por óxido de etileno	Αποστείρωση με αιθυλοξείδιο	تعقيم من خلال أكسيد الإيثيلين	经环氧乙烷灭菌	エチレンオキシド滅菌	Sterilizace ethylenoxidem
	Nicht erneut sterilisieren	Do not sterilize	Ne pas restériliser	Non sterilizzare nuovamente	No volver a esterilizar	Não esterilizar	Δεν επαναπο-τερίωνεται	لا تعقم مرة أخرى	不可重复灭菌	滅菌回しない	Nesterilizujte opakovaně
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	Non riutilizzabile	No reutilizar	Não reutilizar	Δεν επανα-προσποσθεται	غير مسموح بإعادة استخدامها	一次性使用	再使用禁止	Není určeno k opakovanému použití
	Vorsicht! Bitte Warnhinweise beachten	Caution! Please observe warnings!	Attention! Respecter les mises en garde	Attenzione! Si prega di rispettare le avvertenze	¡Tenere en cuenta las instrucciones de uso!	Cuidado! Leia os avisos!	Προσοχή! Λαμβάνεται υπόψη τις υποδείξεις προειδοποίησης	تنبيه! يرجى مراعاة إرشادات الاستخدام	注意！請遵守所有警告說明	注意！全ての警告をお守りください	Pozor! Dodržujte prosím výstražné pokyny
	Gebrauchs-anweisung ist zu beachten!	Operating instructions must be observed!	Respecter le mode d'emploi !	Rispettare le istruzioni per l'uso!	¡Tener en cuenta las instrucciones de uso!	As instruções de operação devem ser respeitadas!	Πρέπει να λαμβάνονται υπόψη οι οδηγίες χρήσης!	يجب مراعاة إرشادات الاستخدام!	務必遵守使用說明！	使用説明書の内容をお守りください！	Dodržte návod k použití!
	Produkt nicht bei Personen anwenden die einen Herzschritt-macher oder implantierten Defibrillator tragen!	Do not use this product on patients with a pacemaker or implanted defibrillator!	Ne pas utiliser le produit sur des personnes portant un stimulateur cardiaque ou un défibrillateur implanté!	Non utilizzare il prodotto per i pazienti portatori di pacemaker o defibrillatore impiantato!	¡El producto no debe ser utilizado por personas que lleven implantado un marcapasso o un desfibrilador!	Não utilize este produto em doentes com marca-passo ou desfibrilador implantado!	Μην χρησιμοποιείτε το προϊόν σε άτομα που φέρουν βηματοδότη ή εμφυτευμένο σπινθηρίστη!	لا يُستخدم المنتج من قبل الأشخاص الذين يعانون من مشاكل في نظم ضربات القلب الطبيعي أو مزود الرجفان القلبي للرج!	不要將產品用於配戴心臟起搏器或植入式除顫器使用方人員。	心臓ペースメーカーまたは植え込み型除細動器を使用されている方の近くで本製品をご使用にならないでください。	Nepoužívejte v blízkosti osob, které nosí kardiostimulátor nebo implantovaný defibrilátor!
	Keine MRT-Untersuchung bei Trägern dieses Produktes!	No MRT examinations on patients using this product!	Ne pas soumettre les personnes portant ce produit à un examen RM !	Non sottoporre i pazienti portatori di questo prodotto ad un esame TRM!	¡No realizar ningún análisis con tomografía por resonancia magnética (TRM) si consume este producto!	Proibido tomografia de ressonância magnética em doentes usando este produto!	Απαγορεύεται κάθε είδους μαγνητικής τομογραφίας για άτομα που φέρουν βηματοδότη ή εμφυτευμένο σπινθηρίστη!	لا يخضع تلقى هذا المنتج إلى فحص التصوير بالرنين المغناطيسي!	不要對配戴該產品的人進行MRI檢查！	本製品を使用されている方の近くでMRI検査を行わないでください！	Neobtěžte tohoto výrobku nesmí vystařeny MRI!
	Bei beschädigter Verpackung nicht verwenden	Do not use if package is damaged	Ne pas utiliser si l'emballage est endommagé	Non usare se la confezione è danneggiata	No utilizar en caso de embalaje dañado	Não utilize se a embalagem estiver danificada	Να μην χρησιμοποιείται εάν η συσκευασία είναι ελαττωματική	لا تستخدم في حالة التلف	包裝有損時不可使用	包裝時有損壞、使用しない下さい	Nepoužívejte, je-li obal poškozený
	Trocken aufbewahren	Keep dry	Conservar au sec	Conservare in luogo asciutto	Conservar en un lugar seco	Mantner seco	Φυλάσσετε σε στεγνό μέρος	يحفظ في مكان جاف	儲存于干燥處	乾燥した場所	Skladujte v suchu
	Oben	Up	Haut	Sopra	Arriba	Em cima	Επάνω	لأعلى	向上	上方	Nahoru
	Temperaturbegrenzung	Temperature limit	Limite de température	Limite di temperatura	Limite de temperatura	Limite de temperatura	Όριο θερμοκρασίας	تحديد درجة الحرارة	溫度限制	溫度制限	Omezení teploty
	Von Sonnenlicht fernhalten	Keep away from sunlight	Tenir éloigné de la lumière du soleil	Tenersi lontano dai raggi del sole	Mantener apartado de la luz solar	Proteger da luz solar	Προσταθείτε από την ηλιακή ακτινοβολία	يُحفظ بعيداً عن ضوء الشمس	避免阳光直射	日光から遠ざけてください	Udržujte z dosahu slunečního svitu
	Zulässige Stapellast	Stacking limit by mass	Charge de gerbage autorisée	Carico di impilaggio consentito	Carga aptable permitida	Carga máxima sobre a embalagem	Επιτρεπτό φορτίο συσκευής	يُسمح بالتحميل عليه	允許最大堆疊質量	重量による積上げ制限	Povolené sňovací zatížení
	Verschreibungs-pflichtig (US-Symbol)	Available only on prescription (US symbol)	Prescription médicale obligatoire (symbole utilisé aux États-Unis)	Soggetto a prescrizione medica (simbolo US)	De prescripción obligatoria (símbolo US)	Disponível somente mediante receita médica (símbolo US)	Υποχρεωτικά συνταγοποιούμενο (συμβολο ΝΗΑ)	يُصرف (رمز الولايات المتحدة)	仅适用于处方治疗 (美国标识)	処方箋 (米 国シンボル)	Pouze na předpis (symbol USA)

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Please note: This medical product may be obtained and used only by trained medical personnel.

1. Description
Pusher and pusher with ring marker

2. Contents and packaging
- Pusher or
- Pusher with ring marker:
- CH 5 = gray color
- CH 6 = black color
- CH 7 = yellow color
- CH 8 = orange color

- Double packaging
- Instructions for use

3. Intended use
Pushers are used to position ureteral stents precisely at the location of implantation via a guide wire.

4. Indication
The medical indication corresponds to that of the ureteral stents.

5. Contraindications
Contraindications result from restrictions for the corresponding basic applications.

6. Possible side effects
There are no known side effects from the use of pushers.

7. Instructions
a) Push the pusher behind the ureteral stent via the guide wire.
b) Only withdraw the guide wire
i. after placing the ureteral stent
ii. then remove the pusher

8. Possible complications and/or risks
- When a hydrophilic wire is used, there is a risk that the pusher will remain stuck to the wire if the wire is no longer adequately hydrated.
- Pushing the stent forward too hard with the pusher can damage the bladder wall.

9. Instructions for use / precautions
- This product is only sterile if the packaging is undamaged and unopened.
- For single use only!
- Do not use stents after the date of expiry has passed.
- Do not wipe the product and set components with alcohol/disinfectant agents. This may damage the surface.
- All components must be carefully checked for compatibility and integrity before use.
- When inserting via the working channel of the cystoscope, the size of the pusher must be adjusted to the working channel.
- Do not damage the pusher with sharp surgical instruments - this increases the risk of breakage. Do not excessively stretch the stent before use. Stretching it too far increases the risk of breakage.
- Pushers are intended for temporary use ≤ 60 minutes.

10. Interaction with other drugs
The pusher could potentially discolor due to the concomitant administration of medication, use of dye-based urine presentation methods, or even as a result of contact with urine.

11. Transport and storage conditions
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

12. Disposal
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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