DE

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EN

## Instructions for use Transurethral Balloon Catheter

IFU Nr 24-A Ballonkatheter EN\_07/21.06.2017

Please note: These medical products may only be obtained and used by trained medical personnel.

- 1. Description**  
Two-way and three-way transurethral balloon catheters made from medical grade silicone or latex with silicone elastomer coating.
- 2. Contents and packaging**  
– Two-way or three-way catheter  
– Catheter stopper (only for 2-way catheters)  
– Instructions for use
- 3. Intended use**  
Temporary urinary diversion, tamponade and possibility of bladder irrigation with 3-way balloon catheters.
- 4. Indications**  
– Diagnostics of the lower urinary tract  
– Intra-operative and post-operative urine drainage  
– Bladder emptying disorders  
– For transurethral bladder catheterization  
– For tamponade and/or lavage following resections in the lower urinary tract

- 5. Contraindications**  
The contraindications known for the respective procedure must be observed. The balloon catheter may not be used in the case of known allergies to the respective catheter material (silicone or latex).

- 6. Possible side effects**  
– Pressure ulceration or inflammation of the urethral mucosa  
– Infections (bacteriuria)  
– Catheter incrustation and obstruction  
– Urethral injury  
– Hematuria

- 7. Instructions**  
When using this product, the general hygiene rules applicable to catheterization must be complied with.  
**Insertion technique (man)**  
The size of the balloon catheter is selected on an individual patient basis in accordance with medical standards. Conventional medical lubricant is drizzled onto the urethral opening using a disposable syringe. The outer urethral meatus is then widened manually, the cone positioned and the lubricant instilled by extending the head. To ensure that the lubricant does not leak out, the urethra should be held compressed for a moment or closed off with a penis clamp. The penis is then held lateral to the glans penis. The catheter is then inserted with the sterile sheathing or sterile forceps. Its free end is grasped between the pinkie and ring finger and secured in this way. On catheters with curved tips, the tip must point upwards. The catheter is then advanced smoothly past the sphincter with constant checks until it is inside the bladder and urine is draining out. The catheter is then advanced a further 5 cm and, using the valve and the Luer attachment on a syringe, blocked through the second lumen and adjusted in the bladder by pulling it backwards.

- Insertion technique (women)**  
The size of the balloon catheter is selected in accordance with medical standards. The labia majora are disinfected using one swab each from the pubic bone towards the anus and then spread using the thumb and index finger. The labia minora are disinfected in the same manner. Throughout the entire catheterization procedure, one hand spreads the labia while catheterization is carried out using only the other hand. Once the lubricant has been instilled, the catheter is introduced into the urethral opening and slowly advanced until it reaches the bladder and urine drains out. The catheter is then advanced a further 5 cm and, using the valve and the Luer attachment on a syringe, blocked through the second lumen and adjusted in the bladder by pulling it backwards. Once the catheter has been inserted, it is connected via the valve to a one-way connection on a urine bag and firmly pressed.

- Insertion technique (child)**  
When selecting the catheter it must be ensured that the catheter is not too tight a fit in the urethra!

- Removal of the catheter**  
– The balloon must be completely drained before removal. To do this, place an empty syringe on the valve and aspirate the fluid. The catheter is then carefully removed.

- 8. Possible complications and/or risks**  
– Long-term catheterization can cause shrinkage of the bladder.  
– Infections  
– Incrustation  
– Catheter dislocation  
– Lack of proper care when advancing the catheter may damage the bladder wall.

- 9. Instructions for use / precautions**  
– This product is only sterile if the packaging is undamaged and unopened.  
– For single use only! Do not re-use, recycle, or re-sterilize. Re-using, recycling, or re-sterilizing can impair the structural integrity of the instrument and/or lead to malfunctioning of the instrument, which can in turn lead to illness, injury, or death of the patient.  
– Do not use catheter after the date of expiry has passed.  
– Do not treat the product with alcohol/disinfectant agents/lubricants. This may damage the surface.  
– All components must be carefully checked for compatibility and integrity before use.  
– Do not over-block the balloon! Only block to the maximum volume specified! (Details can be found on the catheter attachment and/or on the label)  
– Regularly check the patency and position of the catheter. The fill volume of the balloon should be checked at intervals of 1-2 weeks.  
– To reduce the risk of urinary tract infections, the connection between the catheter and the urine bag should not be released. Closed drainage systems are recommended.  
– To refill the balloon, drain it completely first. Fill the balloon with slight pressure. To reduce the risk of loss of volume in the balloon, distilled water + 10% glycerin should be used.  
– In exceptional cases, it may be impossible to unblock the balloon. In this case, the balloon can be emptied through special measures (e.g. by cutting off the valve, piercing the balloon by suprapubic incision, etc.).  
– The user must decide on the most appropriate measure based on the individual situation. Following one of these measures, the balloon must be checked for completeness.  
– Do not clamp the catheter; use a sealing plug if necessary.  
– Silicone balloon catheters must be changed after no more than 30 days.  
– Check that the patient does not have an allergy to latex before using catheters that contain latex. The duration of use is limited to 4 days, and repeated use to 30 days. If long-term urinary catheterization is required, silicone balloon catheters should be used.

- 10. Interaction with other drugs**  
It is possible that the stent may discolor due to the concomitant use of medications or even as a result of contact with urine.

- 11. Transport and storage conditions**  
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

- 12. Disposal**  
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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