IFU Nr. 12-A

Symbols see part B

Instructions for use Pyeloplastic catheter

IFU Nr 12-A Pyeloplastik EN 01/25.02.2016

Please note: This medical product may be obtained and used only by trained medical personnel.

- Description**
Pyeloplastic catheter with pigtail in the center.
- Contents and packaging**
 - Pyeloplastic catheter (biocompatible material, aliphatic PUR)
 - Needle adapter with funnel
 - Double packaging
 - Patient record and instructions for use
- Intended use**
Postoperative urinary drainage after a pyeloplastic operation
- Indications**
 - Congestion in the renal pelvis and ureter after tumor or stone formation
 - Idiopathic sub-pelvic stenosis (intrinsic)
 - High ureteral outflow (intrinsic)
 - Connective tissue symphysis between ureter and lower pole of the kidney (extrinsic)
 - Aberrant vessels (extrinsic)
- Contraindications**
 - Inoperable or weakly functioning kidneys
 - Severe urinary tract infection
 - Hypocoagulability
 - Coagulation disorders
- Possible side effects**
 - Dysuria
 - Flank pain
 - Hematuria
 - Foreign body irritation in the kidney
 - Erosions and ulcerations of the ureter
- Instructions**
Anesthesia and incision of the patient (access to kidney and ureter) are not included in these instructions.
The course of action for the catheter insertion is described in the following.
The pyeloplastic catheter is inserted intra-operatively if the dorsal wall of the new anastomosis is sutured between the ureter and renal pelvis.
For the anastomosis, use a suture material with a longer resorption time (around 35 days)!
 - The ureter side segment (=short end) of the the catheter is placed in the ureter in such a way (control the positioning with graduation!) that the anastomosis bridges the gap securely and the loop portion of the catheter comes to lie tension-free in the new (reduced) renal pelvis.
 - The renal side discharge segment can be fixed on the surface of the kidney using an overcast suture with short-term absorbable suture material (7 - 12 days), in order to avoid dislocation.
 - The pyeloplastic catheter can be implanted through the operation cut or through a special incision out of the abdominal wall.
 - Then the cannula adapter with a funnel can be placed on the catheter and an appropriate drainage system can be connected.
 - This way the function of the kidney can be monitored for several days after the operation. Flushing should only occur in exceptional cases to check the patency (NaCl solution 0.9%)!

* **Alternative:** Another option for the drainage placement is insertion from the renal pelvis. In this case, the bulb-headed probe is attached to the long (externally conveyed) catheter end and is led under visual guidance through the appropriate calyx of the renal pelvis through the parenchyma. When the fibrous capsule is reached, an incision is made. The catheter is pulled into the renal pelvis far enough that the loop segment of the catheter lies tension-free in the renal pelvis.

- Possible complications and/or risks**
 - Dislocation of the catheter.
 - Fragmentation of the catheter with subsequent surgical or endoscopic removal of the product.
 - Encrustation and blockage of the pyeloplastic catheter, rendering it non-functional and potentially causing urinary calculi to form in the kidneys.
 - Post-operative bleeding from the renal parenchyma in the renal pelvis can occur as a result of the transrenal implementation of the drainage system.
 - Strictures can develop due to ureteral damage.
 - Urine loss due to leakage at the puncture side in the renal pelvis.
- Warnings / Precautions**
 - This product is only sterile if the packaging is undamaged and unopened.
 - For single use only!
 - Do not use catheter after the date of expiry has passed.
 - Do not wipe the product and set components with alcohol/disinfectant agents. This may damage the surface.
 - All components must be carefully checked for compatibility and integrity before use.
 - The positioning of the catheter through the renal parenchyma outwards must be done in the least traumatic way possible, resulting in minimal bleeding.
 - Do not damage the catheter with sharp surgical instruments - this increases the risk of breakage.
 - If the patient has a urinary tract infection, the product can become encrusted and blocked relatively quickly.
 - The product is intended for a service life of ≤ 30 days.
 - Check the catheter regularly for drainage effectiveness and position. In order to ensure flawless functioning of the pyeloplastic catheter, regular examinations (urinary status, position of the catheter, etc.) must be performed.
 - The catheter should be laid "tension-free" in order to prevent post-operative bleeding.
 - Please note: Continuous urinary drainage from the drainage system should be particularly ensured if there is only one kidney!
 - If a blockage of the drainage holes (e.g. due to blood clots) is determined post-operatively, it should be flushed out immediately!
 - It must be ensured that the catheter is not also stitched during the suturing of the renal pelvis closure!
- Interaction with other drugs**
The catheter could potentially discolour due to the concomitant administration of medication or even as a result of contact with urine.
- Transport and storage conditions**
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.
- Disposal**
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines of the respective country.

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