IFU Nr. 32-A

EN
Instructions for use
MAGNETIC BLACK STAR / MAGNETIC BLACK STAR-XT

IFU Nr. 32-A Magnetic Black Star EN 04/09.06.2017

Please note: This medical product may only be obtained and used by specialist physicians, physicians and trained personnel.

- Description**
 - Ureteral stent (double-J / single-J) with magnet and accessories
 - Retrieval device (magnetic stent removal system)
- Contents and packaging**
 - Ureteral stent set**
 - Pusher and ureteral stent with bladder-side magnet
 - Guide wire (in sets with guide wire)
 - Patient ID and instructions for use
 - Retrieval Device**
 - Disposable catheter with magnetic Tiemann tip
 - Instructions for use
- Intended use**

Ureteral stents ensure a urinary diversion from the kidney to the bladder, prevent extravasation of urine, or splint the ureter.
- Indications**
 - Stenting of the ureter pre- and postoperatively
 - Prophylactic or postoperative urine drainage security after URS
 - Counteracting of strictures, also scar slits after ureter slits or slits at the ostium
 - Ureteral stenoses of various etiology
 - Ureteral strictures, stenoses, kinks, or adhesion
 - Postoperative, e.g. after Ureteropelvic surgery or ureterotomy
- Contraindications**
 - Serious infections
 - Prolonged obstruction
 - Contracted or spastic bladder
 - Urinary bladder fistulas
 - Urosepsis
 - Urge incontinence
 - Chronic urinary retention or extrinsic obstruction.
 - Cardiac pacemaker
- Possible side effects**
 - Foreign body irritation in the bladder and renal regions
 - Erosions and ulcerations of the ureter
 - Dysuria
 - Hematuria
 - Elevated micturition frequency, nighttime urge to urinate
 - Feeling of pressure in the flank, flank pain
 - Vesicoureteral reflux
 - Incrustation
 - Ascending infections
- Instructions**
 - Pass the guide wire, flexible tip first, into the renal pelvis.
 - Inserting the double-J ureteral stent with magnet**
 - Thread the ureteral stent with the magnet and pusher onto the guide wire. Make sure that the magnet (bladder-side end) is also led over the guide wire.
 - Push the ureteral stent into the renal pelvis via the guide wire using the pusher with x-ray guidance.
 - Inserting the single-J ureteral stent with magnet**
 - Thread the ureteral stent and the pusher, in which the the magnet is pre-mounted through a sleeve, onto the guide wire.
 - Remove the sleeve from the pusher!
 - Push the ureteral stent into the renal pelvis via the guide wire using the pusher with x-ray guidance. While doing so take care that the ostium is only intubated by the tapered part of the shaft of the pusher and that the magnet remains in the bladder.
 - Removing the guide wire / pusher**
 - As soon as the desired position of the ureteral stent has been achieved in the renal pelvis, the guide wire can carefully be withdrawn (the pigtails of the ureteral stent will hereby be formed). Hold the ureteral stent in position using the pusher while the guide wire is being withdrawn.
 - After the successful removal of the guide wire, carefully withdraw the pusher.
 - Removing the ureteral stent with the retrieval device**
 - To remove the ureteral stent, carefully insert the retrieval device into the bladder retrogradely (a retrieval device of size "L" must be used to remove a FR6 or FR7 ureteral stent).
 - Create the connection to the magnet on the ureteral stent by rotating the Tiemann tip (in front of the respective ostium). Successful connection between the retrieval device and the stent magnet is indicated with an audible click.
 - After both magnets are connected, remove the ureteral stent by careful, steady pulling.

- Possible complications and/or risks**
 - Dislocation of the ureteral stent
 - Fragmentation of the stent with subsequent surgical or endoscopic removal of the product
 - Incrustation and blockage of the ureteral stent, rendering it non-functional and potentially causing urinary calculi to form in the kidney
 - Arterial ureteral fistula in patients with vascular reconstruction and stents for chronic vascular disease
 - Strictures of the urethra and ureter
 - Impaired x-ray visibility as a result of air in the lower abdomen of the patient
 - Prevention of guide wire insertion due to obstructions in the ureter
 - Injury to the ureter as a result of improper manipulations of the guide wire
 - Injury to the bladder wall as a result of the application of too much pressure when advancing the ureteral stent with the pusher
 - Injury to the urethra as a result of improper insertion of the magnet (magnet not led using the guide wire)
 - Injury to the urinary tract as a result of the improper positioning of the magnet during insertion (bladder-side end of the ureteral stent in front)
 - Magnetic removal not possible. In this case standard cystoscopic removal must be performed.
- Warnings / Precautions**
 - This product is only sterile if the packaging is undamaged or unopened!
 - For single use only! Do not re-use, recycle, or sterilize multiple times. Re-using, recycling or re-sterilizing the instrument can impair its structural integrity and cause malfunctions that may lead to illness, injury or death of the patient.
 - Do not use stents after the date of expiry has passed.
 - Do not wipe the product and set components with alcohol/disinfectant agents. This may damage the surface
 - All components must be carefully checked for compatibility and integrity before use.
 - Do not damage the stent, pusher, guide wire, and retrieval device with sharp surgical instruments or excessively extend the ureteral stent - this increases the risk of breakage.
 - Carry out a percutaneous nephrostomy where there are extensive defects in the ureter. There is a risk of lesion enlargement from retrograde advancement of the ureteral stent.
 - To ensure correct functioning, a ureteral stent must be selected that corresponds with the length of the ureter.
 - To avoid injury to the bladder wall, ostium, ureter and parenchyma, advance the guide wire, pusher and ureteral stent carefully and in a controlled manner using the magnet.
 - To prevent injury to the urethra, take care that the magnet is threaded onto the guide wire. When using single-J ureteral stents, also take care that the mounting sleeve is removed after threading is complete!
 - If obstructions are detected while inserting the wire, do not increase pressure as this may result in injury.
 - To prevent injuries in the ureteral area, the guide wire must not be led in or out through the pre-positioned ureteral stent (this can cause the wire tip to exit through the drainage holes).
 - The correct placement and complete formation of the pigtails must be ensured using radiological or contrast agent-based fluoroscopy techniques in order to minimize the risk of ureteral stent migration.
 - This product is not a permanent implant.
 - In order to prevent loss of function due to encrustation, the max. service life of a magnetic ureteral stent is limited to 30 days.
 - To detect infections, dislocation or closures early, it is recommended to regularly check urine status and the position and patency of the ureteral stent. In case of irregularities, it is necessary to exchange the ureteral stent.
 - In the event of magnet migration into the ureter, the ureteral stent must be removed cystoscopically.
 - MRT examinations are not permitted for patients with ureteral stents.
 - Patients with cardiac pacemakers may not have the product implanted due to the magnet.
- Interaction with other drugs**

The ureteral stent could potentially discolour due to the concomitant administration of medication, use of dye-based urine presentation methods, or even as a result of contact with urine.
- Transport and storage conditions**

The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5° - 30°C.
- Disposal**

After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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