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# Instructions for Use GUIDE WIRE

GANL Nr. 16\_EN /1/08.04.2015

Please note: This medical product may only be obtained and used by trained medical personnel.

## 1. Description/ Contents

Guide wire are composed of stainless steel with PTFE jacket applied over the core to provide endoscopic visualization of wire movement.

## 2. Indications

Guide wire are designed to be used during urological endoscopic procedures for catheter or other device introduction and exchange.

## 3. Notes for Use and precautionary measures

- This guide wire should be used only by physicians experienced in the use of such devices.
- Sterile unless intimate packaging is damaged or opened
- For single use only; Do not resterilize!
- Sterilized with ethylene oxide.

## 4. Contraindications

Same as for the basic procedure.

## 5. Side effects

Operate guide wire as the procedure of position ureteral stent, otherwise will cause ureteral mucosa damage.

## 6. Anleitung

### Stents one end closed

- Push the stent catheter with the open end in front onto the guidewire (use rigid end!).
- Stent must be fixed on guidewire under tension with clamp.
- Insert stent into the cystoscope channel. When the clamp reaches the opening of the insertion-channel, remove the clamp.
- Push catheter by aid of albarran under x-ray-control into the ureter and the renal pelvis
- When placement is complete, remove guidewire and place the stent by aid of the pusher (blue) in the bladder. Then carefully remove pusher from cystoscope.

### Stents both ends open

- The guidewire with flexible tip in front, should be passed by the obstruction right to the renal pelvis.
- Push the stent catheter over the guidewire with the help of the pusher (blue) and under x-ray control, into the renal pelvis.
- When the stent catheter distal end reaches the renal pelvis stop the advancement. The guidewire has to be removed while the Pusher fixes the stent in its position. The stent will form spontaneously. Then remove the pusher carefully from the cystoscope.

### Intra-operative stents

- Open the ureter at the level of the renal pelvis and push the catheter with the guidewire into the bladder.
- Remove the guidewire and then insert the wire through the drainage hole, which is placed in the middle of the catheter.
- Push the cranial part of the catheter into the renal pelvis.
- After removing the guidewire and checking the correct position of the catheter, the ureter is going to be sewed.

## 7. Interaction with other drugs

Unknown.

## 8. Disposal

After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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