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Instructions for use
Transurethral Balloon Catheter

GANL Nr. 15-A_EN /7/10.09.2019

Please note: This medical product may only be obtained and used by trained nursing and medical personnel. Non-intended use of the product can cause serious and irreversible damage.

- 1. Description**
Two-lumen transurethral balloon catheters made from medical grade silicone or latex with silicone elastomer coating.
- 2. Contents and packaging**
- 2-way standard, Nelaton and Tiemann tip
- 2-way Nelaton tip with integrated balloon
- 2-way Nelaton tip, supersoft
- 2-way latex with silicone elastomer coating, Nelaton tip
- 2-way whistle tip catheter
- Catheter plug and, if necessary, stiffening mandrin
- Double packaging and instructions for use

3. Intended use
Longterm (repeated) transurethral use for draining urine from the bladder.

- 4. Indication**
- Diagnostics of the lower urinary tract
- Bladder emptying disorders
- For transurethral bladder catheterization

- 5. Contraindications**
- Acute prostatitis
- Suspected urethral trauma
- Known allergy to the catheter material

- 6. Possible side effects**
- Pressure ulceration or inflammation of the urethral mucosa
- Infections (bacteriuria)
- Catheter incrustation and obstruction
- Urethral injury
- Hematuria

7. Instructions
When using this product, the general hygiene rules applicable to catheterization must be complied with
Insertion technique (men)
The size of the balloon catheter is selected on an individual patient basis in accordance with medical standards.
Conventional medical lubricant is drizzled onto the urethral opening using a disposable syringe. The outer urethral meatus is then widened manually, the cone positioned and the lubricant instilled by extending the head. To ensure that the lubricant does not leak out, the urethra should be held compressed for a moment or closed off with a penis clamp.
With the left hand, the penis is then held lateral to the glans penis. The catheter is then inserted with the sterile sheathing or sterile forceps. Its free end is grasped between the pinkie and ring finger and secured in this way. On catheters with curved tips, the tip must point upwards. The catheter is then advanced smoothly past the sphincter with constant checks until it is inside the bladder and urine is draining out. The catheter is then advanced by a few cm and, using the valve and the Luer attachment on a syringe, blocked through the second lumen and adjusted in the bladder by pulling it backwards. Once the catheter has been inserted, it is connected via the valve to a one-way connection on a urine bag and firmly pressed.

Insertion technique (women)
The size of the balloon catheter for women is also selected in accordance with medical standards.
The labia majora are disinfected using one swab each from the pubic bone towards the anus and then spread using the thumb and index finger. The labia minora are disinfected in the same manner.
Throughout the entire catheterization procedure, one hand spreads the labia while catheterization is carried out using only the other hand. Once the lubricant has been instilled, the catheter is introduced into the urethral opening and slowly advanced until it reaches the bladder and urine drains out.

The catheter is then advanced by a few cm and, using the valve and the Luer attachment on a syringe, blocked through the second lumen and adjusted in the bladder by pulling it backwards. Once the catheter has been inserted, it is connected via the valve to a one-way connection on a urine bag and firmly pressed.

Removal of the catheter
The balloon must be completely drained before removal. To do this, place an empty syringe on the valve and aspirate the fluid until a vacuum forms. The catheter is then carefully removed.

- 8. Possible complications and/or risks**
- Long-term catheterization can cause shrinkage of the bladder.
- Infections
- Incrustation
- Catheter dislocation
- Lack of proper care when advancing the catheter may damage the bladder wall.
- Bladder stones can damage the balloon.

9. Instructions for use / Precautions
- This product is only sterile if the packaging is undamaged and unopened.
- For single use only! Do not re-use, recycle, or re-sterilize. Re-using, recycling, or re-sterilizing can impair the structural integrity of the instrument and/or lead to malfunctioning of the instrument, which can in turn lead to illness, injury, or death of the patient.
- Do not use catheter after the date of expiry has passed
- Do not treat the product with alcohol/disinfectant agents/lubricants. This may damage the surface.
- All components must be carefully checked before use. Products that may be damaged must not be used.
- Do not over-block the balloon! Only block to the maximum volume specified! (Details can be found on the catheter attachment and/or on the label)
- Check the catheter regularly for drainage effectiveness and position. The fill volume of the balloon should be checked at intervals of 1-2 weeks. If necessary, catheter replacement should be considered.
- To reduce the risk of urinary tract infections, the connection between the catheter and the urine bag should ideally not be released. Closed drainage systems are recommended.
- To refill the balloon, drain it completely first. Fill the balloon with slight pressure. To reduce the risk of loss of volume in the balloon, distilled water/+10% glycerin should be used.
- Fully silicone catheters, silicone catheters with a Visioslide coating and silicone elastomer-coated catheters must be changed after no more than 30 days.
- In exceptional cases, it may be impossible to unblock the balloon. In this situation, various measures can be used to drain the balloon (e.g. cutting the valve, puncturing the balloon with a suprapubic incision, etc.). The user must decide on the most appropriate measure based on the individual situation.
- Do not clamp the catheter; use a sealing plug if necessary.
- For catheters with stiffening mandrins, lubricant shall be applied to the balloon catheter before the mandrin is inserted in order to make removal of the mandrin easier.
- Check that the patient does not have an allergy to latex before using catheters that contain latex.

10. Interaction with other drugs
The catheter could potentially discolor due to the concomitant administration of medication or even as a result of contact with urine.

11. Transport and storage conditions
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport. Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.

12. Disposal
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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