DE

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EN Instructions for use  
Steerable and non-steerable ureteral stent sets

GANL Nr. 01-A Ureterschienen EN\_07/03.04.2018

Please note: These medical products may only be obtained and used by trained medical personnel.

- Description**  
Double-J ureteral stents and stent sets for urinary diversion.
- Contents and packaging**  
- Ureteral stent or  
- tumor stent (with metal spiral or from a hard shaft with soft pigtails)  
- Guide wire GW (PTFE-coated)  
- Pusher or steering device (for steerable sets)  
- Naval clamp (only for ureteral stents with a closed tip)  
- Double packaging  
- Patient record and instructions for use  
- In the steerable sets (Type I and Visiomatic), the stent and steering device are pre-assembled
- Intended use**  
Ureteral stents should ensure a urinary diversion from the kidney to the bladder, prevent extravasation of urine, or splint the ureter.
- Indications**  
- Stenting of the ureter pre- and postoperatively  
- Urinary drainage, even during pregnancy  
- Prophylactic or postoperative urinary drainage security after URS  
- Counteracting of strictures, also scar silts after ureter stents or silts at the ostium  
- Ureteral stenoses of various etiology  
- Ureteral strictures, stenoses, kinks, or adhesion  
- Compression of the ureter due to tumors and other processes which require space in the vicinity  
- Postoperative, e.g. after renal pelvic surgery or ureterotomy
- Contraindications**  
- Serious infections  
- Prolonged obstruction  
- Contracted or spastic bladder  
- Bladder fistulas  
- Dysuria  
- Urge incontinence  
- Chronic urinary retention or extrinsic obstruction.
- Possible side effects**  
- Foreign body irritation in the bladder and renal regions  
- Erosions and ulcerations of the ureter  
- Dysuria  
- Hematuria  
- Elevated micturition frequency, nighttime urge to urinate  
- Feeling of pressure in the flank, flank pain  
- Vesicoureteral reflux  
- Incrustation  
- Ascending infections
- Instructions**  
**Non-steerable sets**  
Ureteral stent renal-side closed.  
a) Insert the guide wire, rigid end first, via the ureteral stent open on the bladder side.  
b) Secure the ureteral stent under tension onto the guide wire using the clamp.  
c) Insert the ureteral stent into the instrument channel of the cystoscope. After reaching the working channel, the naval clamp must be removed.  
d) The ostium is probed and the ureteral stent advanced under fluoroscopic guidance into the renal pelvis.  
e) Following correct placement in the renal pelvis, withdraw the guide wire and lower the ureteral stent into the bladder using the pusher. Then carefully remove the pusher from the cystoscope.  
**Ureteral stent open on both sides**  
a) Advance the guide wire with the flexible tip forward past the obstruction into the renal pelvis.  
b) Under fluoroscopic guidance, push the ureteral stent via the guide wire into the renal pelvis using the pusher.  
c) When the distal end of the ureteral stent reaches the renal pelvis, do not push the ureteral stent any further. The guide wire is withdrawn while the pusher holds the ureteral stent in position. Once the guide wire has been removed, the pigtails of the ureteral stent form.

Typ 1	Visiomatic	Dirigenta
<b>1. Inserting the guide wire (only for open stents)</b> Push the GW forward into the renal pelvis. Avoid bending or loop formation in the bladder when doing so. Wet hydrophilic GW beforehand. For wires with a movable core, the core of the GW can be pulled back, thus making the tip flexible.		
<b>2. Connecting the stent and steering device</b> a) <b>open stents</b> Steering device and ureteral stents are already pre-connected. Thread both together onto the indwelling GW. <b>Note: Do not release the ureteral stent from the inserted GW and connect steering device before use, as it will not be possible to reconnect it.</b> b) <b>closed stents</b> Push the GW forward until the pigtail is stretched out. If the GW is correctly positioned, attach it to the steering device (clamp coupling). When using a wire with movable core, the flexibility of the tip can be individually adjusted (with sizes FR 6 and 7).		

- 3. Insertion into the urinary tract**  
Slowly insert the connected system through the ostium into the ureter via the instrument channel and position it in the renal pelvis.
- 4. Correct positioning of the ureteral stent**  
Pull back the guide wire far enough that the pigtail unrolls and move forward or backward and rotate the ureteral stent to position as desired. Rotational movements and axial movements can be applied to the catheter.
- 5. Releasing the fixing and removing the guide wire**  
After correct placement of the stent, release the clamp coupling of the Y-adapter or the fixing on the bracket and remove the GW.
- 6. Disconnection**  
Slide the outer tube against the inner tube. Hold the handle (red) and twist the tip (blue) to the left as indicated by the arrows. Remove the safeguard and push the bracket together.

- Possible complications and/or risks**  
- Dislocation of the ureteral stent,  
- Fragmentation of the product with subsequent surgical or endoscopic removal  
- Incrustation and blockage of the ureteral stent, rendering it non-functional and potentially causing urinary calculi to form in the kidney,  
- Arterial ureter fistulas in patients with vascular reconstruction and stents in cases of chronic vascular disease  
- Pushing the ureteral stent forward with the guide wire can damage the ureter.  
- Retracting the guide wire and pushing it forward again can cause the tip of the guide wire to exit through the holes in the ureteral stent, causing damage to the ureter area.  
- Pushing the stent forward too hard with the pusher can damage the bladder wall.  
- Due to obstructions in the ureter, it is not possible to advance the guide wire into the kidney,  
- Strictures can develop due to damage to the ureter.  
- Manipulations with the guide wire can trigger ureteral spasms,  
- Damage to the urinary tract when pushing the steering device forward if the stent has already been removed  
- Unintentional premature removal or self-loosening of the ureteral stent
- Warnings / Precautions**  
- This product is only sterile if the packaging is undamaged and unopened.  
- For single use only! Do not re-use, recycle, or re-sterilize. Re-using, recycling, or re-sterilizing can impair the structural integrity of the instrument and/or lead to malfunctioning of the instrument, which can in turn lead to illness, injury, or death of the patient.  
- Do not use the product after the date of expiry has passed.  
- Do not use the product and set components with alcohol/disinfectant agents. This may damage the surface.  
- Check all components for compatibility and integrity before use.  
- When inserting the stent via the working channel of the cystoscope, the size of the ureteral stent must be adjusted to the working channel.  
- Do not damage the set components with sharp surgical instruments – this increases the risk of breakage. Do not excessively stretch the stent before use. Stretching it too far increases the risk of breakage.  
- Correct positioning and complete forming of the pigtail is required to minimize the risk of stent migration. This should be ensured through radiological or contrast agent-assisted fluoroscopy techniques.  
- Any air in the lower abdomen of the patient could impair visibility of the ureteral stent in an x-ray.  
- The ureteral stents are approved for long-term use (> 30 days). However, these are not permanent implants, and it is therefore recommended that the urine status, stent position and patency of the ureteral stent be checked by radiography, isotopic methods or cystoscopy at regular intervals, and at least every 3-4 weeks, to ensure early detection of any infection, dislocation, or blockages which might necessitate a stent change.  
- Where there are extensive defects in the ureter, a percutaneous nephrostomy should be carried out, since blind pushing of the stent incurs the risk of further enlarging any lesion.  
- If a suture is used at the end on the bladder side, this should be removed after a few days. The maximum duration of use for the suture should not exceed 30 days. It should be ensured before usage that the suture can be easily moved.  
- **MRT examinations are not permitted for patients implanted with ureteral stents with metal coils.**  
- **Patients with cardiac pacemakers may not have the product implanted due to the metal coil.**  
- Push the guide wire forward carefully and in a controlled manner to avoid damaging the ostium, the ureter, or the parenchyma.  
- The movable core of the GW must only be advanced in a straight position.  
- The movable core of the guide wire must not be twisted or pushed forward too forcefully, as otherwise the wire coil could be punctured.  
- Do not continue to push the GW forward after withdrawal.  
- Do not continue to push the steering device forward after disconnection.  
- Only remove the safeguard if the ureteral stent has been correctly positioned.  
- Do not use the product if it is clearly damaged.
- Interact with other drugs**  
The stent could potentially dislodge due to the concomitant administration of medication or even as a result of contact with urine.
- Transport and storage conditions**  
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport.  
Products must be stored dry and protected from direct sunlight, in the temperature range of 5-30 degrees Celsius.
- Disposal**  
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines of the respective country.

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