



If you want to have everything under control,
make sure of it:

UROBITE

Single-use Forceps



**Special Design
for reliable grip**



Our disposable UROBITE forceps in reusable quality combines the design of alligator mouth and rat tooth. UROBITE for the safe removal of stents and foreign bodies is distinguished by their ease of use, functionality and efficiency.

- Practical disposable instrument that saves on pre-treatment, sterilisation and validation costs
- Ergonomic grip for comfort and control during the procedure
- Two five-toothed rows and rat tooth for reliable gripping
- Highly flexible tube for protecting the endoscope
- Compatible with all popular cystoscopes
- Minimise the risk of cross contamination

Specialist for
endourological products

No compromise on stent removal

The UROBITE single-use forceps expands our range of high quality endourology medical devices. We value usability. That means: When removing ureteral splints or foreign bodies, the surgeon always works with a new, guaranteed sterile instrument, whose features can support the procedure optimally. Decide on our cost-effective UROBITE forceps and reduce time requirement, costs and risks of reprocessing without having to concern about decreasing material quality.

Our benchmark for working comfortably and cost-effectively:

- Single-use forceps in reusable quality
- Highly flexible tube in two sizes for application in rigid and flexible cystoscopes
- Can be used with all common cystoscopes with a working channel of at least 6 FR
- Stainless steel forceps with powerful five-toothed rows and rat tooth
- Comfort and control through the ergonomic, easy-to-use handle
- Time and efficiency advantage by eliminating the reprocessing process



Ordering Information:

Order number	Diameter mm / FR	Length	Shape	SU
GSF-18-600	1.8 mm	60 cm	Alligator forceps with rat tooth	10
GSF-18-1200	1.8 mm	120 cm	Alligator forceps with rat tooth	10

We are urology specialists and your reliable partner when it comes to optimal user and patient comfort. Ask us, and we will offer you the best product for your needs.

Distributed by:

UROTECH GmbH

Medi-Globe-Straße 1-5
D-83101 Rohrdorf OT Achenmühle
Phone: +49 (0) 8032 973-210
Fax: +49 (0) 8032 973-211
E-Mail: info@urotech.com
www.urotech.com

Manufacturer:
Anrei Medical (HZ) Co., Ltd.