


# We are URO



Specialist for  
endourological  
products

## CONTENTS



URETERAL STENTS  
Chapter 1



NEPHROSTOMY  
Chapter 2



STONE BASKETS  
Chapter 3




URETERAL CATHETERS  
Chapter 4



SURGICAL CATHETERS  
Chapter 5



URODYNAMICS  
Chapter 6



WIRES  
Chapter 7



ACCESSORIES  
Chapter 8



BIOPSY  
Chapter 9



SUPRAPUBIC BLADDER  
DRAINAGE  
Chapter 10



TRANSURETHRAL  
CATHETERS  
Chapter 11



URINE DRAINAGE  
SYSTEMS  
Chapter 12



**UROVISION<sup>®</sup>**  
A MEDI-GLOBE COMPANY


**URETERAL STENTS**

Ureteral Stents, standard	1.1
Ureteral Stents, soft	1.6
Tumor Stents	1.10
ESWL Ureteral Stent Set	1.14
VISIOSTAR Intraoperative	1.15
VISIOSTENT Stenosis	1.16
VISIOFLUX Single-J Ureteral Stent set	1.17
Intubation Stent	1.18
Pyeloplasty Catheters	1.19
VISIO WEB	1.20
VISIO-WEB SOFT	1.21
Complete Sets	1.22

## URETERAL STENTS

# VISIOSTAR STANDARD Double-J Ureteral Stent Set

### DESCRIPTION



#### Ureteral Stent

- ◆ Polyurethane
- ◆ Drainage openings in pigtails and shaft
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

#### Pusher

- ◆ Length: 40 cm

#### Guide wire

- ◆ PTFE coating

#### Clip

- ◆ In sets with closed tip



## URETERAL STENTS

### VISIOSTAR STANDARD Double-J Ureteral Stent Set

Open tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-194722	ST-1954722	4.7	22	0.028"	150	straight	standard
ST-194724	ST-1954724	4.7	24	0.028"	150	straight	standard
ST-194726	ST-1954726	4.7	26	0.028"	150	straight	standard
ST-194728	ST-1954728	4.7	28	0.028"	150	straight	standard
ST-194730	ST-1954730	4.7	30	0.028"	150	straight	standard
ST-190624	ST-195624	6	24	0.035"	150	straight	standard
ST-190626	ST-195626	6	26	0.035"	150	straight	standard
ST-190628	ST-195628	6	28	0.035"	150	straight	standard
ST-190630	ST-195630	6	30	0.035"	150	straight	standard
ST-190724	ST-195724	7	24	0.035"	150	straight	standard
ST-190726	ST-195726	7	26	0.035"	150	straight	standard
ST-190728	ST-195728	7	28	0.035"	150	straight	standard
ST-190730	ST-195730	7	30	0.035"	150	straight	standard
ST-190824	ST-195824	8	24	0.038"	150	straight	standard
ST-190826	ST-195826	8	26	0.038"	150	straight	standard
ST-190828	ST-195828	8	28	0.038"	150	straight	standard
ST-190830	ST-195830	8	30	0.038"	150	straight	standard
ST-190924	ST-195924	9	24	0.038"	150	straight	standard
ST-190926	ST-195926	9	26	0.038"	150	straight	standard
ST-190928	ST-195928	9	28	0.038"	150	straight	standard
ST-190930	ST-195930	9	30	0.038"	150	straight	standard
ST-191024	ST-1951024	10	24	0.045"	150	straight	standard
ST-191026	ST-1951026	10	26	0.045"	150	straight	standard
ST-191028	ST-1951028	10	28	0.045"	150	straight	standard
ST-191030	ST-191030	10	30	0.045"	150	straight	standard

Variants with sutures on the bladder end are available on request.

## URETERAL STENTS

Closed tip						
Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-184722	4.7	22	0.028"	150	straight	standard
ST-184724	4.7	24	0.028"	150	straight	standard
ST-184726	4.7	26	0.028"	150	straight	standard
ST-184728	4.7	28	0.028"	150	straight	standard
ST-184730	4.7	30	0.028"	150	straight	standard
ST-180624	6	24	0.035"	150	straight	standard
ST-180626	6	26	0.035"	150	straight	standard
ST-180628	6	28	0.035"	150	straight	standard
ST-180630	6	30	0.035"	150	straight	standard
ST-180724	7	24	0.035"	150	straight	standard
ST-180726	7	26	0.035"	150	straight	standard
ST-180728	7	28	0.035"	150	straight	standard
ST-180730	7	30	0.035"	150	straight	standard
ST-180824	8	24	0.038"	150	straight	standard
ST-180826	8	26	0.038"	150	straight	standard
ST-180828	8	28	0.038"	150	straight	standard
ST-180830	8	30	0.038"	150	straight	standard
ST-180924	9	24	0.038"	150	straight	standard
ST-180926	9	26	0.038"	150	straight	standard
ST-180928	9	28	0.038"	150	straight	standard
ST-180930	9	30	0.038"	150	straight	standard
ST-181024	10	24	0.045"	150	straight	standard
ST-181026	10	26	0.045"	150	straight	standard
ST-181028	10	28	0.045"	150	straight	standard
ST-181030	10	30	0.045"	150	straight	standard


Variants with sutures on the bladder end are available on request.

## URETERAL STENTS

# VISIOSTAR STANDARD Double-J Ureteral Stent Set

Flexible length

### DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Multiple coils on both sides for shaft length flexibility from 24 to 32 cm
- ◆ Drainage openings in pigtails and shaft
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque



### Pusher

- ◆ Length: 40 cm

### Guide wire

- ◆ PTFE coating

### Clip

- ◆ In sets with closed tip





## URETERAL STENTS

### Open tip

Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-394724-32	ST-3974724-32	4.7	24-32	0.028"	150	straight	standard
ST-390624-32	ST-397624-32	6	24-32	0.035"	150	straight	standard
ST-390724-32	ST-397724-32	7	24-32	0.035"	150	straight	standard
ST-390824-32	ST-397824-32	8	24-32	0.038"	150	straight	standard


### Closed tip

Order number with Guide wire		Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-384724-32		4.7	24-32	0.028"	150	straight	standard
ST-380624-32		6	24-32	0.035"	150	straight	standard
ST-380724-32		7	24-32	0.035"	150	straight	standard
ST-380824-32		8	24-32	0.038"	150	straight	standard

## URETERAL STENTS

# VISIOSTAR STANDARD Double-J Ureteral Stent Set Steerable

### DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Drainage openings in pigtails and shaft
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Guide unit

- ◆ Pusher
- ◆ Disconnecter
- ◆ Length: 53 cm

### Guide wire

- ◆ PTFE coating

### Clip

- ◆ In sets with closed tip





## URETERAL STENTS


Open tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire (inch)	Wire length (cm)	Wire tip	Wire stiffness
US-110524	US-120524	4.7	24	0.028"	150	straight	standard
US-110526	US-120526	4.7	26	0.028"	150	straight	standard
US-110528	US-120528	4.7	28	0.028"	150	straight	standard
US-110530	US-120530	4.7	30	0.028"	150	straight	standard
US-110624	US-120624	6	24	0.035"	150	straight	standard/movable core
US-110626	US-120626	6	26	0.035"	150	straight	standard/movable core
US-110628	US-120628	6	28	0.035"	150	straight	standard/movable core
US-110630	US-120630	6	30	0.035"	150	straight	standard/movable core
US-110724	US-120724	7	24	0.035"	150	straight	standard/movable core
US-110726	US-120726	7	26	0.035"	150	straight	standard/movable core
US-110728	US-120728	7	28	0.035"	150	straight	standard/movable core
US-110730	US-120730	7	30	0.035"	150	straight	standard/movable core
US-110824	US-120824	8	24	0.038"	150	straight	standard
US-110826	US-120826	8	26	0.038"	150	straight	standard
US-110828	US-120828	8	28	0.038"	150	straight	standard
US-110830	US-120830	8	30	0.038"	150	straight	standard
US-110924	US-120924	9	24	0.038"	150	straight	standard
US-110926	US-120926	9	26	0.038"	150	straight	standard
US-110928	US-120928	9	28	0.038"	150	straight	standard
US-110930	US-120930	9	30	0.038"	150	straight	standard

Closed tip							
Order number with Guide wire		Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
US-100524		4.7	24	0.028"	150	straight	standard
US-100526		4.7	26	0.028"	150	straight	standard
US-100528		4.7	28	0.028"	150	straight	standard
US-100530		4.7	30	0.028"	150	straight	standard
US-100624		6	24	0.035"	150	J-tip	standard/movable core
US-100626		6	26	0.035"	150	J-tip	standard/movable core
US-100628		6	28	0.035"	150	J-tip	standard/movable core
US-100630		6	30	0.035"	150	J-tip	standard/movable core
US-100724		7	24	0.035"	150	J-tip	standard/movable core
US-100726		7	26	0.035"	150	J-tip	standard/movable core
US-100728		7	28	0.035"	150	J-tip	standard/movable core
US-100730		7	30	0.035"	150	J-tip	standard/movable core

## URETERAL STENTS

# VISIOMATIC® Double-J Ureteral Stent Set Steerable with VISIOMATIC® Discharge System

### DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Drainage openings in pigtails and shaft
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Visiomatic discharge system

- ◆ Pusher
- ◆ Disconnecter
- ◆ Y-adapter with crimp coupling for securing the Guide wire
- ◆ Length: 51 cm

### Guide wire

- ◆ PTFE coating




## URETERAL STENTS

Open tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
US-110524-V	US-120524-V	4.7	24	0.028"	150	straight	standard
US-110526-V	US-120526-V	4.7	26	0.028"	150	straight	standard
US-110528-V	US-120528-V	4.7	28	0.028"	150	straight	standard
US-110530-V	US-120530-V	4.7	30	0.028"	150	straight	standard
US-110624-V	US-120624-V	6	24	0.035"	150	straight	standard/movable core
US-110626-V	US-120626-V	6	26	0.035"	150	straight	standard/movable core
US-110628-V	US-120628-V	6	28	0.035"	150	straight	standard/movable core
US-110630-V	US-120630-V	6	30	0.035"	150	straight	standard/movable core
US-110724-V	US-120724-V	7	24	0.035"	150	straight	standard/movable core
US-110726-V	US-120726-V	7	26	0.035"	150	straight	standard/movable core
US-110728-V	US-120728-V	7	28	0.035"	150	straight	standard/movable core
US-110730-V	US-120730-V	7	30	0.035"	150	straight	standard/movable core
US-110824-V	US-120824-V	8	24	0.038"	150	straight	standard
US-110826-V	US-120826-V	8	26	0.038"	150	straight	standard
US-110828-V	US-120828-V	8	28	0.038"	150	straight	standard
US-110830-V	US-120830-V	8	30	0.038"	150	straight	standard
US-110924-V	US-120924-V	9	24	0.038"	150	straight	standard
US-110926-V	US-120926-V	9	26	0.038"	150	straight	standard
US-110928-V	US-120928-V	9	28	0.038"	150	straight	standard
US-110930-V	US-120930-V	9	30	0.038"	150	straight	standard

## URETERAL STENTS

# MONOSTAR STANDARD Mono-J Ureteral Stent Set

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Drainage openings only available for pigtail
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Guide wire

- ◆ PTFE coating

### Funnel adapter

- ◆ For CH/FR 04 to 09

### Clip

- ◆ In sets with closed tip





## URETERAL STENTS

### Mono-J with Double Pigtail

Order number with Guide wire		Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-250700-DP		7	70	0.035"	150	straight	standard
ST-250800-DP		8	70	0.038"	150	straight	standard
ST-250900-DP		9	70	0.038"	150	straight	standard

### Open tip

Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-250400	ST-257400	4	70	0.025"	150	straight	standard
ST-250500	ST-257500	5	70	0.028"	150	straight	standard
ST-250600	ST-257600	6	70	0.035"	150	straight	standard
ST-250700	ST-257700	7	70	0.035"	150	straight	standard
ST-250800	ST-257800	8	70	0.038"	150	straight	standard
ST-250900	ST-257900	9	70	0.038"	150	straight	standard
ST-251000		10	70	0.045"	150	straight	standard
ST-251200		12	70	0.045"	150	straight	standard


### Closed tip

Order number with Guide wire		Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-240470		4	70	0.025"	150	straight	standard
ST-240570		5	70	0.028"	150	straight	standard
ST-240670		6	70	0.035"	150	straight	standard
ST-240770		7	70	0.035"	150	straight	standard
ST-240870		8	70	0.038"	150	straight	standard

## SOFT URETERAL STENTS

### VISIOSTAR SOFT Double-J Ureteral Stent Set Soft material

#### DESCRIPTION



#### Ureteral Stent

- ◆ Polyurethane
- ◆ Especially soft material to ensure a high level of patient comfort
- ◆ Drainage openings in pigtails and shaft
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

#### Pusher

- ◆ Length: 40 cm

#### Guide wire

- ◆ PTFE coating

#### Clip

- ◆ In sets with closed tip







Open tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-214722	ST-2174722	4.7	22	0.028"	150	straight	standard
ST-214724	ST-2174724	4.7	24	0.028"	150	straight	standard
ST-214726	ST-2174726	4.7	26	0.028"	150	straight	standard
ST-214728	ST-2174728	4.7	28	0.028"	150	straight	standard
ST-214730	ST-2174730	4.7	30	0.028"	150	straight	standard
ST-210624	ST-217624	6	24	0.035"	150	straight	standard
ST-210626	ST-217626	6	26	0.035"	150	straight	standard
ST-210628	ST-217628	6	28	0.035"	150	straight	standard
ST-210630	ST-217630	6	30	0.035"	150	straight	standard
ST-210724	ST-217724	7	24	0.035"	150	straight	standard
ST-210726	ST-217726	7	26	0.035"	150	straight	standard
ST-210728	ST-217728	7	28	0.035"	150	straight	standard
ST-210730	ST-217730	7	30	0.035"	150	straight	standard
ST-210824	ST-217824	8	24	0.038"	150	straight	standard
ST-210826	ST-217826	8	26	0.038"	150	straight	standard
ST-210828	ST-217828	8	28	0.038"	150	straight	standard
ST-210830	ST-217830	8	30	0.038"	150	straight	standard


Closed tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-204722	ST-2054722	4.7	22	0.028"	150	straight	standard
ST-204724	ST-2054724	4.7	24	0.028"	150	straight	standard
ST-204726	ST-2054726	4.7	26	0.028"	150	straight	standard
ST-204728	ST-2054728	4.7	28	0.028"	150	straight	standard
ST-204730	ST-2054730	4.7	30	0.028"	150	straight	standard
ST-200624	ST-205624	6	24	0.035"	150	straight	standard
ST-200626	ST-205626	6	26	0.035"	150	straight	standard
ST-200628	ST-205628	6	28	0.035"	150	straight	standard
ST-200630	ST-205630	6	30	0.035"	150	straight	standard
ST-200724	ST-205724	7	24	0.035"	150	straight	standard
ST-200726	ST-205726	7	26	0.035"	150	straight	standard
ST-200728	ST-205728	7	28	0.035"	150	straight	standard
ST-200730	ST-205730	7	30	0.035"	150	straight	standard
ST-200824	ST-205824	8	24	0.038"	150	straight	standard
ST-200826	ST-205826	8	26	0.038"	150	straight	standard
ST-200828	ST-205828	8	28	0.038"	150	straight	standard
ST-200830	ST-205830	8	30	0.038"	150	straight	standard

Variants with sutures on the bladder end are available on request.

## SOFT URETERAL STENTS

### VISIOSTAR SOFT Double-J Ureteral Stent Set Steerable, soft material

#### DESCRIPTION



#### Ureteral Stent

- ◆ Polyurethane
- ◆ Especially soft material to ensure a high level of wearing comfort
- ◆ Drainage openings in pigtails and shaft
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

#### Guide unit

- ◆ Pusher
- ◆ Disconnecter
- ◆ Length: 53 cm

#### Guide wire

- ◆ PTFE coating

#### Clip

- ◆ In sets with closed tip





## SOFT URETERAL STENTS


Open tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
US-140524	US-150524	4.7	24	0.028"	150	straight	standard
US-140526	US-150526	4.7	26	0.028"	150	straight	standard
US-140528	US-150528	4.7	28	0.028"	150	straight	standard
US-140530	US-150530	4.7	30	0.028"	150	straight	standard
US-140624	US-150624	6	24	0.035"	150	straight	standard/movable core
US-140626	US-150626	6	26	0.035"	150	straight	standard/movable core
US-140628	US-150628	6	28	0.035"	150	straight	standard/movable core
US-140630	US-150630	6	30	0.035"	150	straight	standard/movable core
US-140724	US-150724	7	24	0.035"	150	straight	standard/movable core
US-140726	US-150726	7	26	0.035"	150	straight	standard/movable core
US-140728	US-150728	7	28	0.035"	150	straight	standard/movable core
US-140730	US-150730	7	30	0.035"	150	straight	standard/movable core
US-140824	US-150824	8	24	0.038"	150	straight	standard
US-140826	US-150826	8	26	0.038"	150	straight	standard
US-140828	US-150828	8	28	0.038"	150	straight	standard
US-140830	US-150830	8	30	0.038"	150	straight	standard

## SOFT URETERAL STENTS

# VISIOMATIC® Double-J Ureteral Stent Set

Soft material, steerable with VISIOMATIC® discharge system

### DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Especially soft material to ensure a high level of patient comfort
- ◆ Drainage openings in pigtails and shaft
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Visiomatic discharge system

- ◆ Pusher
- ◆ Disconnecter
- ◆ Y-adaptor with crimp coupling for securing the Guide wire
- ◆ Length: 51 cm

### Guide wire

- ◆ PTFE coating



## SOFT URETERAL STENTS


Open tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
US-140524-V	US-150524-V	4.7	24	0.028"	150	straight	standard
US-140526-V	US-150526-V	4.7	26	0.028"	150	straight	standard
US-140528-V	US-150528-V	4.7	28	0.028"	150	straight	standard
US-140530-V	US-150530-V	4.7	30	0.028"	150	straight	standard
US-140624-V	US-150624-V	6	24	0.035"	150	straight	standard/movable core
US-140626-V	US-150626-V	6	26	0.035"	150	straight	standard/movable core
US-140628-V	US-150628-V	6	28	0.035"	150	straight	standard/movable core
US-140630-V	US-150630-V	6	30	0.035"	150	straight	standard/movable core
US-140724-V	US-150724-V	7	24	0.035"	150	straight	standard/movable core
US-140726-V	US-150726-V	7	26	0.035"	150	straight	standard/movable core
US-140728-V	US-150728-V	7	28	0.035"	150	straight	standard/movable core
US-140730-V	US-150730-V	7	30	0.035"	150	straight	standard/movable core
US-140824-V	US-150824-V	8	24	0.038"	150	straight	standard
US-140826-V	US-150826-V	8	26	0.038"	150	straight	standard
US-140828-V	US-150828-V	8	28	0.038"	150	straight	standard
US-140830-V	US-150830-V	8	30	0.038"	150	straight	standard

## SOFT URETERAL STENTS

### MONOSTAR SOFT Mono-J Ureteral Stent Set

Soft material

#### DESCRIPTION



#### Ureteral Stent

- ◆ Polyurethane
- ◆ Especially soft material to ensure a high level of patient comfort
- ◆ Drainage openings only available for pigtail
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

#### Guide wire

- ◆ PTFE coating

#### Adapter


- ◆ For attaching to the urine drainage bag

Open tip							
Order number with Guide wire	Order number without Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-350600	ST-357600	6	70	0.035"	150	straight	standard
ST-350700	ST-357700	7	70	0.035"	150	straight	standard
ST-350800	ST-357800	8	70	0.038"	150	straight	standard
ST-350900	ST-357900	9	70	0.038"	150	straight	standard

## SPECIAL STENTS

# SPECIAL VISIOSTENT TUMOR STENT Double-J Ureteral Stent Set

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Produced from a single piece, so it is not possible for the pigtail to loosen
- ◆ Stable shaft for maximal pressure resistance
- ◆ Soft pigtails for improved wearing comfort
- ◆ No drainage openings in the shaft, to prevent infiltration of tumor tissue growth
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Pusher

- ◆ Length: 40 cm

### Guide wire

- ◆ PTFE coating

## SPECIAL STENTS

# SPECIAL VISIOSTENT TUMOR STENT Double-J Ureteral Stent Set

Open tip						
Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire $\sigma$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
TS-290624	6	24	0.035"	150	straight	standard
TS-290626	6	26	0.035"	150	straight	standard
TS-290628	6	28	0.035"	150	straight	standard
TS-290630	6	30	0.035"	150	straight	standard
TS-290724	7	24	0.035"	150	straight	standard
TS-290726	7	26	0.035"	150	straight	standard
TS-290728	7	28	0.035"	150	straight	standard
TS-290730	7	30	0.035"	150	straight	standard
TS-290824	8	24	0.038"	150	straight	standard
TS-290826	8	26	0.038"	150	straight	standard
TS-290828	8	28	0.038"	150	straight	standard
TS-290830	8	30	0.038"	150	straight	standard
TS-290924	9	24	0.038"	150	straight	standard
TS-290926	9	26	0.038"	150	straight	standard
TS-290928	9	28	0.038"	150	straight	standard
TS-290930	9	30	0.038"	150	straight	standard


Variants with sutures on the bladder end are available on request.



## SPECIAL STENTS

# VISIOMATIC<sup>®</sup> TUMOR STENT Double-J Ureteral Stent Set Steerable with VISIOMATIC<sup>®</sup> Discharge System

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Produced from a single element
- ◆ Stable shaft for maximal pressure resistance
- ◆ Soft pigtails for improved wearing comfort
- ◆ No drainage openings in the shaft, to prevent infiltration of tumor tissue growth
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Visiomatic Discharge System

- ◆ Pusher
- ◆ Disconnecter
- ◆ Y-adapter with crimp coupling for securing the Guide wire
- ◆ Length: 51 cm

### Guide wire

- ◆ PTFE coating

## SPECIAL STENTS


# VISIOMATIC TUMOR STENT Double-J Ureteral Stent Set Steerable with VISIOMATIC® Discharge System

Open tip						
Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
TS-110624-V	6	24	0.035"	150	straight	standard
TS-110626-V	6	26	0.035"	150	straight	standard
TS-110628-V	6	28	0.035"	150	straight	standard
TS-110630-V	6	30	0.035"	150	straight	standard
TS-110724-V	7	24	0.035"	150	straight	standard
TS-110726-V	7	26	0.035"	150	straight	standard
TS-110728-V	7	28	0.035"	150	straight	standard
TS-110730-V	7	30	0.035"	150	straight	standard
TS-110824-V	8	24	0.038"	150	straight	standard
TS-110826-V	8	26	0.038"	150	straight	standard
TS-110828-V	8	28	0.038"	150	straight	standard
TS-110830-V	8	30	0.038"	150	straight	standard
TS-110924-V	9	24	0.038"	150	straight	standard
TS-110926-V	9	26	0.038"	150	straight	standard
TS-110928-V	9	28	0.038"	150	straight	standard
TS-110930-V	9	30	0.038"	150	straight	standard

## SPECIAL STENTS

# ONKOVISION<sup>®</sup> TUMOR STENT Double-J Ureteral Stent Set With a stainless steel helix

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Highly flexible stainless steel helix for maximal pressure resistance with minimal ureter irritation
- ◆ Soft pigtails for improved wearing comfort
- ◆ No drainage openings in the shaft, to prevent infiltration of tumor tissue growth
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Pusher

- ◆ Length: 40 cm

### Guide wire

- ◆ PTFE coating

## SPECIAL STENTS

# ONKOVISION® TUMOR STENT Double-J Ureteral Stent Set With a stainless steel helix


Open tip						
Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
TS-310624	6	24	0.032"	150	straight	standard
TS-310626	6	26	0.032"	150	straight	standard
TS-310628	6	28	0.032"	150	straight	standard
TS-310630	6	30	0.032"	150	straight	standard
TS-310724	7	24	0.035"	150	straight	standard
TS-310726	7	26	0.035"	150	straight	standard
TS-310728	7	28	0.035"	150	straight	standard
TS-310730	7	30	0.035"	150	straight	standard
TS-310824	7	24	0.035"	150	straight	standard
TS-310826	8	26	0.035"	150	straight	standard
TS-310828	8	28	0.035"	150	straight	standard
TS-310830	8	30	0.035"	150	straight	standard

## SPECIAL STENTS

# ONKOVISION<sup>®</sup> TUMOR STENT Double-J Ureteral Stent Set

With a stainless steel helix, Steerable with VISIOMATIC<sup>®</sup> discharge system

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Highly flexible stainless steel helix in the shaft for maximal pressure resistance with minimal ureter irritation
- ◆ Soft pigtails for improved wearing comfort
- ◆ No drainage openings in the shaft, to prevent infiltration of tumor tissue growth
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Visiomatic Discharge System

- ◆ Pusher
- ◆ Disconnecter
- ◆ Y-adapter with crimp coupling for securing the Guide wire
- ◆ Length: 51 cm

### Guide wire

- ◆ PTFE coating
- ◆ Guide wire: stiff

## SPECIAL STENTS


### ONKOVISION® TUMOR STENT Double-J Ureteral Stent Set With a stainless steel helix, Steerable with VISIOMATIC® discharge system

Open tip						
Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
TS-330624	6	24	0.032"	150	straight	standard
TS-330626	6	26	0.032"	150	straight	standard
TS-330628	6	28	0.032"	150	straight	standard
TS-330630	6	30	0.032"	150	straight	standard
TS-330724	7	24	0.035"	150	straight	standard
TS-330726	7	26	0.035"	150	straight	standard
TS-330728	7	28	0.035"	150	straight	standard
TS-330730	7	30	0.035"	150	straight	standard
TS-330824	8	24	0.035"	150	straight	standard
TS-330826	8	26	0.035"	150	straight	standard
TS-330826	8	28	0.035"	150	straight	standard
TS-330830	8	30	0.035"	150	straight	standard

## SPECIAL STENTS

# VISIOSTAR ESWL Double-J Ureteral Stent Set With special profile

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Special profile to support discharge of stones
- ◆ Drainage openings only in pigtails
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Pusher

- ◆ Length: 40 cm

### Guide wire

- ◆ PTFE coating

### Clip

- ◆ In sets with closed tip



### Open tip

Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-230724	7	24	0.035"	150	straight	standard
ST-230726	7	26	0.035"	150	straight	standard
ST-230728	7	28	0.035"	150	straight	standard
ST-230730	7	30	0.035"	150	straight	standard


### Closed tip

Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-220724	7	24	0.035"	150	straight	standard
ST-220726	7	26	0.035"	150	straight	standard
ST-220728	7	28	0.035"	150	straight	standard
ST-220730	7	30	0.035"	150	straight	standard

## SPECIAL STENTS


# VISIOSTAR INTRAOPERATIVE Double-J Ureteral Stent Set

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Shaft opening central or throughout
- ◆ Closed tip on both ends
- ◆ Radiopaque



### Pusher

- ◆ Length: 40 cm

### Guide wire

- ◆ PTFE coating



## SPECIAL STENTS

### Central opening for Guide wire – closed tip

Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-260524	4.7	24	0.028"	100	straight	standard
ST-260526	4.7	26	0.028"	100	straight	standard
ST-260528	4.7	28	0.028"	100	straight	standard
ST-260530	4.7	30	0.028"	100	straight	standard
ST-260624	6	24	0.035"	100	straight	standard
ST-260626	6	26	0.035"	100	straight	standard
ST-260628	6	28	0.035"	100	straight	standard
ST-260630	6	30	0.035"	100	straight	standard
ST-260724	7	24	0.035"	100	straight	standard
ST-260726	7	26	0.035"	100	straight	standard
ST-260728	7	28	0.035"	100	straight	standard
ST-260730	7	30	0.035"	100	straight	standard
ST-260824	8	24	0.038"	100	straight	standard
ST-260826	8	26	0.038"	100	straight	standard
ST-260828	8	28	0.038"	100	straight	standard
ST-260830	8	30	0.038"	100	straight	standard
ST-260924	9	24	0.038"	100	straight	standard
ST-260926	9	26	0.038"	100	straight	standard
ST-260928	9	28	0.038"	100	straight	standard
ST-260930	9	30	0.038"	100	straight	standard

### Multiple openings for Guide wire – closed tip

Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-270524	4.7	24	0.028"	100	straight	standard
ST-270526	4.7	26	0.028"	100	straight	standard
ST-270528	4.7	28	0.028"	100	straight	standard
ST-270530	4.7	30	0.028"	100	straight	standard
ST-270624	6	24	0.035"	100	straight	standard
ST-270626	6	26	0.035"	100	straight	standard
ST-270628	6	28	0.035"	100	straight	standard
ST-270630	6	30	0.035"	100	straight	standard
ST-270724	7	24	0.035"	100	straight	standard
ST-270726	7	26	0.035"	100	straight	standard
ST-270728	7	28	0.035"	100	straight	standard
ST-270730	7	30	0.035"	100	straight	standard
ST-270824	8	24	0.038"	100	straight	standard
ST-270826	8	26	0.038"	100	straight	standard
ST-270828	8	28	0.038"	100	straight	standard
ST-270830	8	30	0.038"	100	straight	standard
ST-270924	9	24	0.038"	100	straight	standard
ST-270926	9	26	0.038"	100	straight	standard
ST-270928	9	28	0.038"	100	straight	standard
ST-270930	9	30	0.038"	100	straight	standard


Additional sizes available on request.

## SPECIAL STENTS

# VISIOSTENT STENOSIS Double-J Ureteral Stent Set

Partial shaft reinforcement

### DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Partial reinforcement (proximal, central, distal) for dilatation of stenosis
- ◆ Drainage openings only in pigtails
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Pusher

- ◆ Length: 40 cm

### Guide wire

- ◆ PTFE coating



## SPECIAL STENTS

### Bladder-end reinforcement – open tip

Order number with Guide wire	Size (CH/FR)	Shaft reinforcement (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
TP-291226	7	12	26	0.038"	150	straight	standard
TP-291228	7	12	28	0.038"	150	straight	standard
TP-291230	7	12	30	0.038"	150	straight	standard

### Central reinforcement – open tip

Order number with Guide wire	Size (CH/FR)	Shaft reinforcement (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
TM-291226	7	12	26	0.038"	150	straight	standard
TM-291228	7	12	28	0.038"	150	straight	standard
TM-291230	7	12	30	0.038"	150	straight	standard


### Kidney-end reinforcement – open tip

Order number with Guide wire	Size (CH/FR)	Shaft reinforcement (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
TD-291226	7	12	26	0.038"	150	straight	standard
TD-291228	7	12	28	0.038"	150	straight	standard
TD-291230	7	12	30	0.038"	150	straight	standard

## SPECIAL STENTS

### **VISIOFLUX** Single-J Ureteral Stent set For reducing bladder irritation

#### DESCRIPTION



#### Ureteral Stent

- ◆ Polyurethane
- ◆ Thanks to the optimal length of the stent, which ends before the ostium, irritations in the bladder and ostium are prevented
- ◆ Pigtail only at kidney end/proximal
- ◆ With suture at the straight end, which stays in the bladder
- ◆ Drainage openings in shaft and pigtail
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

#### Pusher

- ◆ Length: 40 cm

#### Guide wire

- ◆ PTFE coating

## SPECIAL STENTS

### Open tip

Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
FL-930622	6	22	0.035"	150	straight	standard
FL-930624	6	24	0.035"	150	straight	standard
FL-930626	6	26	0.035"	150	straight	standard
FL-930722	7	22	0.035"	150	straight	standard
FL-930724	7	24	0.035"	150	straight	standard
FL-930726	7	26	0.035"	150	straight	standard
FL-930822	8	22	0.038"	150	straight	standard
FL-930824	8	24	0.038"	150	straight	standard
FL-930826	8	26	0.038"	150	straight	standard

### Closed tip


Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire $\phi$ (inch)	Wire length (cm)	Wire tip	Wire stiffness
FL-920622	6	22	0.035"	150	straight	standard
FL-920624	6	24	0.035"	150	straight	standard
FL-920626	6	26	0.035"	150	straight	standard
FL-920722	7	22	0.035"	150	straight	standard
FL-920724	7	24	0.035"	150	straight	standard
FL-920726	7	26	0.035"	150	straight	standard
FL-920822	8	22	0.038"	150	straight	standard
FL-920824	8	24	0.038"	150	straight	standard
FL-920826	8	26	0.038"	150	straight	standard

## SPECIAL STENTS

# INTUBATION STENT Double-J Ureteral Stent Set

Partial shaft thickening, also in the pigtails

### DESCRIPTION



#### Ureteral Stent

- ◆ Polyurethane
- ◆ Partial reinforcement (proximal, distal) for dilatation of stenosis
- ◆ Drainage openings only in pigtails
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

#### Pusher

- ◆ Length: 40 cm

#### Guide wire

- ◆ PTFE coating

#### Bladder-end reinforcement – open tip – for antegrade application

Order number with Guide wire	Size (CH/FR)	Shaft reinforcement (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
SI-301226	7	12	26	0.038"	150	straight	standard
SI-301228	7	12	28	0.038"	150	straight	standard
SI-301230	7	12	30	0.038"	150	straight	standard

#### Kidney-end reinforcement – open tip – for retrograde application


Order number with Guide wire	Size (CH/FR)	Shaft reinforcement (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
SI-291226	7	12	26	0.038"	150	straight	standard
SI-291228	7	12	28	0.038"	150	straight	standard
SI-291230	7	12	30	0.038"	150	straight	standard

## SPECIAL STENTS

# PYELOPLASTY CATHETERS

For pyeloplasty

## DESCRIPTION



### Catheter

- ◆ Polyurethane
- ◆ Drainage openings in pigtail and ureter section
- ◆ With straight tip or pigtail
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Adapter

- ◆ For attachment to the urine drainage bag

#### With straight tip

Order number	Size (CH/FR)	Length (cm)	Length of ureter section (cm)	Length of drainage section (cm)	Maximum Guide wire ø
UR-910675	6	75	20	50	0.035"
UR-910775	7	75	20	50	0.035"
UR-910875	8	75	20	50	0.038"

#### With pigtail

Order number	Size (CH/FR)	Length (cm)	Length of ureter section (cm)	Length of drainage section (cm)	Maximum Guide wire ø
UR-930628	6	75	28	35	0.035"
UR-930728	7	75	28	35	0.035"

## SPECIAL STENTS

# VISIO-WEB STANDARD

For supraventricular urinary diversion

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Transparent material with x-ray contrast strips and plastic stylet
- ◆ 5 drainage openings at the proximal end for unimpeded urine discharge
- ◆ Open or closed tip
- ◆ Markings for exact positioning
- ◆ Application e.g. following neobladder formation



### Adapter

- ◆ Funnel adapter for rinsing or attachment of urine drainage bag

Closed tip			
Order number	Size (CH/FR)	Length (cm)	Maximum Guide wire ø
WB-106	6	80	0.035"
WB-107	7	80	0.035"
WB-108	8	80	0.038"
WB-109	9	80	0.038"
WB-110	10	80	0.045"

Open tip			
Order number	Size (CH/FR)	Length (cm)	Maximum Guide wire ø
WB-906	6	80	0.035"
WB-907	7	80	0.035"
WB-908	8	80	0.038"
WB-909	9	80	0.045"
WB-910	10	80	0.045"



## SPECIAL STENTS

# VISIO-WEB SOFT

For supraventricular urinary diversion

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Soft material with x-ray contrast strips and plastic stylet
- ◆ 5 drainage openings at the kidney end for unimpeded urine discharge
- ◆ Closed tip
- ◆ Markings for exact positioning
- ◆ Application e.g. following neobladder formation



### Adapter


- ◆ Funnel adapter for rinsing or attachment of urine drainage bag

Closed tip			
Order number	Size (CH/FR)	Length (cm)	Maximum Guide wire ø
WB-206	6	80	0.038"
WB-207	7	80	0.038"
WB-208	8	80	0.038"

## SPECIAL STENTS

# VISIOSTAR COMPLETE SET Double-J Ureteral Stent Set with Ureteral Catheter

## DESCRIPTION



### Ureteral Stent

- ◆ Polyurethane
- ◆ Drainage openings in pigtails and shaft
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque

### Ureteral Catheter

- ◆ Made of Pebax with stainless steel stylet
- ◆ Open Tiemann tip
- ◆ Size: CH/FR 4.7
- ◆ Length: 70 cm

### Pusher

- ◆ Length: 40 cm

### Guide wire

- ◆ PTFE coating

### Adapter

- ◆ Cannula adapter

## SPECIAL STENTS

# VISIOSTAR COMPLETE SET Double-J Ureteral Stent Set with Ureteral Catheter

Open tip						
Order number with Guide wire	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
ST-186524	4.7	24	0.035"	150	straight	standard
ST-186526	4.7	26	0.035"	150	straight	standard
ST-186528	4.7	28	0.035"	150	straight	standard
ST-186530	4.7	30	0.035"	150	straight	standard
ST-186624	6	24	0.035"	150	straight	standard
ST-186626	6	26	0.035"	150	straight	standard
ST-186628	6	28	0.035"	150	straight	standard
ST-186630	6	30	0.035"	150	straight	standard
ST-186724	7	24	0.035"	150	straight	standard
ST-186726	7	26	0.035"	150	straight	standard
ST-186728	7	28	0.035"	150	straight	standard
ST-186730	7	30	0.035"	150	straight	standard

Also available with sutures on the bladder end.

WWW.UROVISION.DE




**NEPHROSTOMY**

RENOVISION Nephrostomy Sets, with Pigtail Catheter	2.1
RENOVISION 2-Step Method	2.2
RENOVISION Nephrostomy, Pigtail Catheters with funnel	2.3
RENOVISION Nephrostomy, Balloon Catheters & Sets	2.4
REINSERTION Catheter Set	2.7
VISIOFIX Renal Cyst Puncture Set	2.8

## NEPHROSTOMY

# RENOVISION NEPHROSTOMY Puncture and Exchange Set With pigtail catheter

## SET COMPONENTS



### Nephrostomy Catheter

- ◆ Polyurethane
- ◆ Drainage openings in pigtail
- ◆ Luer Lock Connection
- ◆ Optimal radiopaque visibility due to sandblasted tip

### Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 18 G, 19.5 cm in length

### Guide wire (Lunderquist)

- ◆ PTFE coating

### Universal Adapter

- ◆ Luer lock connection and silicone tube attachment (polyurethane coating)
- ◆ Length: 10 cm

## NEPHROSTOMY

# RENOVISION NEPHROSTOMY Puncture and Exchange Set with pigtail catheter

Single catheter		
Order number	Size (CH/FR)	Length (cm)
RN-480635	6	30
RN-480645	6	40
RN-480735	7	30
RN-480745	7	40
RN-480835	8	30
RN-480845	8	40
RN-480935	9	30
RN-480945	9	40
RN-481045	10	40
RN-481445	14	40

Puncture Set							
Order number	Size (CH/FR)	Length (cm)	Dilators	Wire ø (inch)	Wire length (cm)	Wire tip (r=1.5 mm)	Wire stiffness
RN-440635	6	30	05/06/07	0.035"	80	J-tip	stiff
RN-440645	6	40	05/06/07	0.035"	80	J-tip	stiff
RN-440735	7	30	06/07/08	0.035"	80	J-tip	stiff
RN-440745	7	40	06/07/08	0.035"	80	J-tip	stiff
RN-440835	8	30	07/08/09	0.035"	80	J-tip	stiff
RN-440845	8	40	07/08/09	0.035"	80	J-tip	stiff
RN-440935	9	30	08/09/10	0.035"	80	J-tip	stiff
RN-440945	9	40	08/09/10	0.035"	80	J-tip	stiff
RN-441045	10	40	08/10/12	0.035"	80	J-tip	stiff

Exchange Set							
Order number	Size (CH/FR)	Length (cm)	Dilators	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
RN-450635	6	30	05/06/07	0.035"	80	straight	stiff
RN-450645	6	40	05/06/07	0.035"	80	straight	stiff
RN-450735	7	30	06/07/08	0.035"	80	straight	stiff
RN-450745	7	40	06/07/08	0.035"	80	straight	stiff
RN-450835	8	30	07/08/09	0.035"	80	straight	stiff
RN-450845	8	40	07/08/09	0.035"	80	straight	stiff
RN-450935	9	30	08/09/10	0.035"	80	straight	stiff
RN-450945	9	40	08/09/10	0.035"	80	straight	stiff
RN-451045	10	40	08/10/12	0.035"	80	straight	stiff




## NEPHROSTOMY

# RENOVISION 2-STEP TECHNIQUE

## With pigtail catheter

### SET COMPONENTS



### Nephrostomy Catheter

- ◆ Polyurethane
- ◆ Drainage openings in pigtail
- ◆ Stretching tool to facilitate insertion of mandrin and puncture needle
- ◆ Luer Lock Connection
- ◆ Optimal radiopaque visibility due to sandblasted tip

### Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 18 G, 19.5 cm in length

### Mandrin

- ◆ Rounded tip to enable safe stretching of catheter

### Guide wire

- ◆ Schüller Guide wire – uncoated stainless steel

### Universal Adapter

- ◆ Luer lock connection and silicone tube attachment (parylene coating)
- ◆ Length: 10 cm


### Two-step method

Order number	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip (r=1.5 mm)	Wire stiffness
RN-900835	8	30	0.035"	80	J-tip	stiff
RN-900935	9	30	0.035"	80	J-tip	stiff

## NEPHROSTOMY

# RENOVISION NEPHROSTOMY Puncture and Exchange Set with Funnel

## SET COMPONENTS



### Nephrostomy Catheter

- ◆ Polyurethane
- ◆ Drainage openings in pigtail
- ◆ Funnel for connecting urine bag
- ◆ Optimal radiopaque visibility due to sandblasted tip

### Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 18 G, 19.5 cm in length

### Guide wire (Lunderquist)

- ◆ PTFE coating

## NEPHROSTOMY

### Single catheter

Order number	Size (CH/FR)	Length (cm)
RN-501145	11	40
RN-501445	14	40

### Puncture Set

Order number	Size (CH/FR)	Length (cm)	Dilators with split shaft*	Wire ø (inch)	Wire length (cm)	Wire tip (r=1.5 mm)	Wire stiffness
RN-461140	11	40	08/10/12*	0.035"	80	J-tip	stiff
RN-461440	14	40	10/12/14/16*	0.035"	80	J-tip	stiff
RN-461640	16	40	10/12/14/16/18*	0.035"	80	J-tip	stiff
RN-461840	18	40	10/14/18/20*	0.035"	80	J-tip	stiff

### Exchange Set

Order number	Size (CH/FR)	Length (cm)	Dilators with split shaft*	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
RN-471140	11	40	08/10/12*	0.035"	80	straight	stiff
RN-471440	14	40	10/12/14/16*	0.035"	80	straight	stiff


## NEPHROSTOMY

# RENOVISION-INTEGRAL NEPHROSTOMY

## Puncture and Exchange Set

Seldinger method, 100 % silicone, with integrated balloon

### SET COMPONENTS



### Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Opening behind the balloon for caliceal drainage
- ◆ Integrated balloon
- ◆ X-ray contrast strips
- ◆ Radiopaque, reinforced tip
- ◆ Markings for exact positioning
- ◆ Optimal radiopaque visibility due to sandblasted tip

### Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 18 G, 19.5 cm in length

### Guide wire (Lunderquist)

- ◆ PTFE coating



## NEPHROSTOMY

### Single catheter

Order number	Size (CH/FR)	Length (cm)	Balloon (ml)
RN-431026	10	30	1.5
RN-421040	10	40	1.5
RN-431226	12	30	2
RN-421240	12	40	2
RN-431426	14	30	3
RN-421440	14	40	3
RN-431626	16	30	5
RN-421640	16	40	5
RN-431826	18	30	5
RN-421840	18	40	5
RN-432026	20	30	5
RN-422040	20	40	5
RN-422240	22	40	5
RN-422440	24	40	5
RN-422640	26	40	5

### Puncture Set

Order number	Size (CH/FR)	Length (cm)	Balloon (ml)	Dilators* with split shaft	Wire (inch)	Wire length (cm)	Wire tip (r=1.5 mm)	Wire stiffness
RN-401040	10	40	1.5	06/08/10/12*	0.035"	80	J-tip	stiff
RN-401240	12	40	2	08/10/12/14*	0.035"	80	J-tip	stiff
RN-401440	14	40	3	10/12/14/16*	0.035"	80	J-tip	stiff
RN-401640	16	40	5	10/12/14/16/18*	0.035"	80	J-tip	stiff
RN-401840	18	40	5	10/14/18/20*	0.035"	80	J-tip	stiff
RN-402040	20	40	5	10/14/18/22*	0.035"	80	J-tip	stiff
RN-402240	22	40	5	10/14/18/22/24*	0.035"	80	J-tip	stiff

### Puncture Set with 3-part Puncture Cannula

Order number	Size (CH/FR)	Length (cm)	Balloon (ml)	Dilators* with split shaft	Wire ø (inch)	Wire length (cm)	Wire tip (r=1.5 mm)	Wire stiffness
RN-401250	12	40	2	08/10/12/14*	0.035"	80	J-tip	stiff
RN-401450	14	40	3	08/10/12/14/16*	0.035"	80	J-tip	stiff

### Exchange Set

Order number	Size (CH/FR)	Length (cm)	Balloon (ml)	Dilators* with split shaft	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
RN-411040	10	40	1.5	06/08/10/12*	0.035"	80	straight	stiff
RN-411240	12	40	2	08/10/12/14*	0.035"	80	straight	stiff
RN-411440	14	40	3	10/12/14/16*	0.035"	80	straight	stiff
RN-411640	16	40	5	12/14/16/18*	0.035"	80	straight	stiff


## NEPHROSTOMY

# RENOVISION-INTEGRAL NEPHROSTOMY

## Puncture Set

Seldinger method, 100 % silicone, with integrated balloon

### SET COMPONENTS



### Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Opening behind the balloon for caliceal drainage
- ◆ Integrated balloon
- ◆ X-ray contrast strips
- ◆ Radiopaque, reinforced tip
- ◆ Markings for exact positioning

### Puncture Needle

- ◆ 2-part needle with trocar tip
- ◆ Medical stainless steel
- ◆ Optimal radiopaque visibility due to sandblasted tip
- ◆ 18 G, 19.5 cm in length

### Guide wire (Lunderquist)

- ◆ PTFE coating


Puncture Set								
Order number	Size (CH/FR)	Length (cm)	Balloon (ml)	Dilators* with split shaft	Wire $\phi$ (inch)	Wire length (cm)	Wire tip (r=1.5 mm)	Wire stiffness
RN-401240-TR	12	40	2	06/08/10/12*	0.035"	80	J-tip	stiff

## NEPHROSTOMY

# RENOVISION BALLOON CATHETER

100 % silicone, with hollow mandrin, integrated balloon

### SET COMPONENTS



### Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Opening behind the balloon for caliceal drainage
- ◆ Integrated balloon
- ◆ Open tip reinforced with metal
- ◆ Markings for exact positioning



### Hollow Mandrin

- ◆ For stabilization of the balloon catheter

### Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 18 G, 19.5 cm in length

### Guide wire (Lunderquist)

- ◆ PTFE coating

## NEPHROSTOMY

# RENOVISION BALLOON CATHETER

100 % silicone, with hollow mandrin, integrated balloon

Single catheter			
Order number	Size (CH/FR)	Length (cm)	Balloon (ml)
RN-511225	12	25	2
RN-511235	12	40	2
RN-511425	14	25	3
RN-511435	14	40	3
RN-511625	16	25	5
RN-511635	16	40	5
RN-511825	18	25	5
RN-512025	20	25	5
RN-512225	22	25	5

Puncture Set								
Order number	Size (CH/FR)	Length (cm)	Balloon (ml)	Dilators	Wire ø (inch)	Wire length (cm)	Wire tip (r=1.5 mm)	Wire stiffness
RN-741235	12	40	2	08/10/12	0.035"	80	J-tip	stiff
RN-741435	14	40	3	08/10/12/14	0.035"	80	J-tip	stiff




## NEPHROSTOMY

# REINSERTION CATHETER SET

For probing of nephrostomy canal after dislocation of nephrostomy catheter

### SET COMPONENTS



#### Catheter

- ◆ made of Pebax
- ◆ Angled tip (45°)
- ◆ Open tip
- ◆ Radiopaque

#### Guide wire (only in set)


- ◆ PTFE coating

Single catheter			
Order number	Size (CH/FR)	Length (cm)	Maximum Guide wire $\phi$
RN-520650	06	45	0.035"

Set with wire						
Order number	Size (CH/FR)	Length (cm)	Wire length (cm)	Wire diameter $\phi$	Wire tip	Wire stiffness
RN-520651	06	45	150	0.035"	straight	standard

# VISIOFIX RENAL CYST PUNCTURE SET with pigtail

## SET COMPONENTS



### Catheter (with pigtail)

- ◆ Polyurethane, fixed
- ◆ Holes in pigtails
- ◆ Open tip
- ◆ Stretching tool to facilitate insertion of puncture needle in the pigtail catheter
- ◆ Optimal radiopaque visibility due to sandblasted tip

### Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 19.5 G

Puncture Set with Pigtail Catheter			
Order number	Size (CH/FR)	Length (cm)	Length of puncture needle (cm)
RN-530525	05	25	26.5



**STONE BASKETS**

## STONE BASKETS

### VISIOCATCH STONE BASKETS

Super-elastic Nitinol core enclosed in a highly elastic polyamide tube\*

#### "TIPLESS"



- ◆ Nitinol core enclosed
- ◆ Flexible tip for good maneuverability in the renal pelvis region
- ◆ High tensile strength, and good pressure and torque stability resulting from lasered baskets
- ◆ Handle can be disconnected and reconnected
- ◆ Optimal kink-resistance due to metal-reinforced inner tube
- ◆ Ergonomic scissor grip for easy handling
- ◆ The basket for the greatest demands

Order number	Size (CH/FR)	Wires	Operating length (cm)	Basket length (mm)	Basket diameter (mm)
DO-99194120	1.9	4	120	15	11

### VISIOCATCH STONE BASKETS

Super-elastic Nitinol

#### "HELICAL"



Order number	Size (CH/FR)	Wires	Operating length (cm)	Basket length (mm)	Basket diameter (mm)	Tip length (mm)
DO-903475	3.0	4	90	40	12	3
DO-9034120	3.0	4	120	40	12	3


## STONE BASKETS

# NiTAN STONE BASKET

## High-quality stone retrieval baskets

Our NiTAN high-quality stone retrieval baskets are designed for capturing and extracting urinary stones from the urinary tract and minimise the risk of injury to patients. Tipless stone baskets are produced using a special laser procedure and are perfect for use with a flexible ureteroscope to extract stones from the kidneys.

### „TIPLESS” AND „HELICAL”





NiTAN helical and tipless stone baskets are available in a range of designs.

- ◆ Ergonomic UV-design handle for a secure hold
- ◆ Super-elastic, excellent kink resistance and flexibility for greater scope
- ◆ High radial force and perfect 1:1 transmission of rotational movements
- ◆ Blue colour coding for excellent endoscopic visibility
- ◆ Supplied in sterile blister pack, ready-for-use

Order Number	Size (CH/FR)	Wires	Working length (cm)	Basket ø (mm)	Basket shape
48012456	2.5	4	90	12	helical
48012457	3.0	4	90	12	helical
48012458	3.0	4	120	12	helical
48012461	4.5	4	90	15	helical
48012462	4.5	5	90	15	helical
48012459	1.8	4	120	11	Tipless
48012460	2.5	4	90	12	Tipless

Commercial products distributed by UROVISION. Not available in all markets!

[WWW.UROVISION.DE](http://WWW.UROVISION.DE)




**URETERAL CATHETER**

## URETERAL CATHETER

### VISION URETERAL CATHETER

For retrograde urography



#### Nelaton tip


Order number Closed tip	Order number Open tip	Size (CH/FR)	Length (cm)	Maximum wire ø (inch)
UR-820370	UR-830370	3	70	0.025"
UR-820470	UR-830470	4	70	0.035"
UR-724770	UR-734770	4.7	70	0.035"
UR-820570	UR-830570	5	70	0.035"
UR-820670	UR-830670	6	70	0.038"
UR-820770	UR-830770	7	70	0.038"

#### Tiemann tip

Order number Closed tip	Order number Open tip	Size (CH/FR)	Length (cm)	Maximum wire ø (inch)
UR-840370	UR-850370	3	70	0.025"
UR-840470	UR-850470	4	70	0.035"
UR-744770	UR-754770	4.7	70	0.035"
UR-840570	UR-850570	5	70	0.035"
UR-840670	UR-850670	6	70	0.038"
UR-840770	UR-850770	7	70	0.038"



## URETERAL CATHETER



### Ureteral catheter by Prof. Walz

- ◆ Two-sided tip configuration
- ◆ Open Nelaton tip
- ◆ Conically opening tip by Prof. Walz
- ◆ Ideal for retrograde imaging of upper urinary tract
- ◆ Optimal sealing of ureter orifice
- ◆ Top visibility in x-ray
- ◆ CM gradation for ideal placement
- ◆ Incl. UK adapter

Ureteral catheter by Prof. Walz				
Order number Closed tip	Size (CH/FR)	Tip by Prof. Walz (CH/ FR)	Length (cm)	Maximum wire ø (inch)
UR-630470	04	07	70	0.035"

Ureter Dilation Set for Ureteral Stenosis				
Order number Closed tip	Sizes in Set (CH/FR)	Length per catheter (cm)	Maximum wire ø (inch)	Extra
HB-5-10	5-10	70	0.035"	Incl. universal adapter and Tuohy Borst adapter

[WWW.UROVISION.DE](http://WWW.UROVISION.DE)



**SURGICAL  
CATHETERS**


VISIODRAIN Balloon Catheters, Silicone	5.1
VISIOCATH Balloon Catheters, Silicone-coated	5.3
Tamponade Catheters	5.5
NEOVISION Neobladder Catheters	5.6

## SURGICAL CATHETERS

# VISIODRAIN BALLOON CATHETER

## Postoperative Tamponade, 100 % silicone

### DESCRIPTION



### Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Large balloon volume for safe compression and effective hemostasis
- ◆ Very high lavage and drainage performance levels
- ◆ 2-way catheter
- ◆ Variable balloon inflation

#### 2-way Balloon Catheters with Couvelaire tip

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-491830-50	18	42	30 - 50
KW-492040-60	20	42	40 - 60
KW-492250-70	22	42	50 - 70
KW-492460-80	24	42	60 - 80

#### 2-way Balloon Catheters with Dufour tip

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-481830-50	18	42	30 - 50
KW-482040-60	20	42	40 - 60
KW-482250-70	22	42	50 - 70
KW-482460-80	24	42	60 - 80

## SURGICAL CATHETERS

# VISIODRAIN BALLOON CATHETER

Postoperative Tamponade, 100 % silicone


2-way Balloon Catheters with Mercier tip			
Order number with open tip	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-521830-50-SL	18	42	30 - 50
KW-522040-60-SL	20	42	40 - 60
KW-522250-70-SL	22	42	50 - 70
KW-522460-80-SL	24	42	60 - 80

## SURGICAL CATHETERS

# VISIODRAIN BALLOON CATHETER

Postoperative irrigation and Tamponade, 100 % silicone

## DESCRIPTION



### Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Large balloon volume for safe compression and effective hemostasis
- ◆ Very high lavage and drainage performance levels
- ◆ 3-way catheter
- ◆ Variable balloon inflation

#### 3-way Balloon Catheters with Nelaton tip

Order number with closed tip	Order number with open tip	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-471630	KW-471630-SL	16	42	30
KW-471830	KW-471830-SL	18	42	30
KW-472030	KW-472030-SL	20	42	30
KW-472230	KW-472230-SL	22	42	30
KW-472430	KW-472430-SL	24	42	30

#### 3-way Balloon Catheters with Couvelaire tip

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-511830-50	18	42	30 - 50
KW-512040-60	20	42	40 - 60
KW-512250-70	22	42	50 - 70
KW-512460-80	24	42	60 - 80

## SURGICAL CATHETERS

# VISIODRAIN BALLOON CATHETER

Postoperative irrigation and Tamponade, 100 % silicone

### 3-way Balloon Catheters with Dufour tip

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-501830-50	18	42	30 - 50
KW-502040-60	20	42	40 - 60
KW-502250-70	22	42	50 - 70
KW-502460-80	24	42	60 - 80

### 3-way Balloon Catheters with Mercier tip

Order number with closed tip	Order number with open tip	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-531830-50	KW-531830-50-SL	18	40	30 - 50
KW-532040-60	KW-532040-60-SL	20	40	40 - 60
KW-532250-70	KW-532250-70-SL	22	40	50 - 70
KW-532460-80	KW-532460-80-SL	24	40	60 - 80




## SURGICAL CATHETERS

# VISIOCATH BALLOON CATHETER

Postoperative tamponade, with latex core and solid silicone coating inside and outside

## DESCRIPTION



## Balloon Catheter

- ◆ Silicone coating for optimal tissue compatibility
- ◆ High balloon stability, which largely eliminates the possibility of balloon ruptures and subsequent blockage
- ◆ Soft balloon catheter offers wearing comfort for patients
- ◆ The catheter is supported by a plastic spiral in the latex core to prevent collapse of the catheter lumen during bladder lavage
- ◆ 2-way catheter

### 2-way Balloon Catheters with Nelaton tip

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-011880	18	40	30 - 50
KI-012080	20	40	40 - 60

### 2-way Balloon Catheters with Couvelaire tip


Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-021880	18	40	30 - 50
KI-022080	20	40	40 - 60
KI-022280	22	40	50 - 70

## SURGICAL CATHETERS

# VISIOCATH BALLOON CATHETER

Postoperative lavage, with latex core and solid silicone coating inside and outside

## DESCRIPTION



## Balloon Catheter

- ◆ Silicone coating for optimal tissue compatibility
- ◆ High balloon stability, which largely eliminates the possibility of balloon ruptures and subsequent blockage
- ◆ Soft balloon catheter offers comfort of use for patients
- ◆ The catheter is supported by a plastic spiral in the latex core to prevent collapse of the catheter lumen during bladder lavage
- ◆ 3-way catheter

### 3-way Balloon Catheters with Nelaton Tip and End Hole

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-061860	18	40	30 - 50
KI-062060	20	40	40 - 60
KI-062260	22	40	50 - 70
KI-062460	24	40	60 - 80

### 3-way Balloon Catheters with Couvelaire tip

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-091880	18	40	30 - 50
KI-092080	20	40	40 - 60
KI-092280	22	40	50 - 70

## SURGICAL CATHETERS

### 3-way Balloon Catheters with Dufour tip

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-071880	18	40	30 - 50
KI-072080	20	40	40 - 60
KI-072280	22	40	50 - 70
KI-072480	24	40	60 - 80

### 3-way Balloon Catheters with Mercier tip


Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-081880	18	40	30 - 50
KI-082080	20	40	40 - 60
KI-082280	22	40	50 - 70
KI-082480	24	40	60 - 80

## SURGICAL CATHETERS

### UROVISION TAMPONADE CATHETERS

With double balloon, with latex core and solid silicone coating inside and outside

#### DESCRIPTION



#### Tamponade catheter

- ◆ Silicone coating for optimal tissue compatibility
- ◆ High balloon stability, which largely eliminates the possibility of balloon ruptures and subsequent blockage
- ◆ Soft balloon catheter offers wearing comfort for patients


Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KD-212240	22	40	1. Balloon 15 ml / 2. Balloon 30 ml

## SURGICAL CATHETERS

# NEOVISION NEOBLADDER CATHETERS

Bladder catheter for postoperative urinary diversion for neobladders following radical prostatectomy

### DESCRIPTION



### Neobladder Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Very low encrustation levels
- ◆ Catheter can be exchanged through the open tip using a wire
- ◆ Catheter does not collapse during strong aspiration of coagulated mucus or blood

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
NB-502040	20	40	30
NB-502240	22	40	30

[WWW.UROVISION.DE](http://WWW.UROVISION.DE)



**URODYNAMIC**

Urodynamic 2-way Measurement Catheters	6.1
Urodynamic 2-way Measurement Catheter with pigtail	6.2
Urodynamic 3-way Measurement Catheters	6.3
Rectal Pressure Measurement Catheters	6.5
Accessories	6.6



## URODYNAMIC

# URODYNAMIC 2-WAY MEASUREMENT CATHETERS

## DESCRIPTION



### 2-way Measurement Catheters

- ◆ PVC-free and latex-free (material: polyurethane)
- ◆ Easy application due to good steering stability
- ◆ Color-coded attachments
- ◆ Radiopaque
- ◆ Optimal measurement results due to special catheter configuration
- ◆ Transurethral and suprapubic bladder pressure measurement

Nelaton tip		
Order number	Size (CH/FR)	Length (cm)
UR-325	07	25
UR-340	07	40
UR-425	09	25
UR-440	09	40


Tiemann tip		
Order number	Size (CH/FR)	Length (cm)
UR-530	07	40
UR-540	09	40

## URODYNAMIC

# URODYNAMIC 2-WAY MEASUREMENT CATHETER

## With pigtail

### DESCRIPTION



### 2-way Measurement Catheter

- ◆ PVC-free and latex-free (material: polyurethane)
- ◆ Easy application due to good steering stability
- ◆ Color-coded attachments
- ◆ Radiopaque
- ◆ Optimal measurement results due to special catheter configuration
- ◆ With stylet and stretching shell
- ◆ Transurethral and suprapubic bladder pressure measurement

With pigtail		
Order number	Size (CH/FR)	Length (cm)
UR-640	09	40

## URODYNAMIC

# URODYNAMIC 3-WAY MEASUREMENT CATHETERS

For measurement of urethra pressure profile

## DESCRIPTION



### 3-way Measurement Catheters

- ◆ PVC-free and latex-free (material: polyurethane)
- ◆ Easy insertion
- ◆ Color-coded attachments
- ◆ Radiopaque
- ◆ Optimal measurement results due to special catheter configuration

Nelaton tip		
Order number	Size (CH/FR)	Length (cm)
UR-140	08	40
UR-240	10	40

Tiemann tip		
Order number	Size (CH/FR)	Length (cm)
UR-150	08	40
UR-250	10	40

## URODYNAMIC

# URODYNAMIC 3-WAY MEASUREMENT CATHETERS

For measurement of urethra pressure profile, with spacers

## DESCRIPTION



### 3-way Measurement Catheters

- ◆ Optimal measurement results through the use of special spacers on the sensor
- ◆ PVC-free and latex-free (material: polyurethane)
- ◆ Easy insertion
- ◆ Color-coded attachments
- ◆ Radiopaque

Nelaton tip		
Order number	Size (CH/FR)	Length (cm)
UR-145	08	40
UR-245	10	40

Tiemann tip		
Order number	Size (CH/FR)	Length (cm)
UR-155	08	40
UR-255	10	40

## URODYNAMIC

# SINGLE-WAY RECTAL PRESSURE MEASUREMENT CATHETERS

Measurement catheter for measuring rectal pressure, 100% silicone

## DESCRIPTION



### 3-way Measurement Catheters


- ◆ Rectal catheter made of 100% silicone
- ◆ Fixing balloon prevents the catheter from slipping out of the rectum in case of compromised pelvic floor muscle tone
- ◆ The end of the catheter is protected from stool by a protective barrier. At the same time, perforations in the protective barrier guarantee a virtually direct connection to the pressure in the rectum once filled with water
- ◆ This ensures maximum precision of measurements
- ◆ The pressure measurement channel can be ventilated up to 100%

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
UR-460	10	30	15
UR-470	18	30	30 - 45

## URODYNAMIC

### ACCESSORIES


#### Three-part extension set



Order number	Inner Diameter (mm)	Length (cm)
UD-Z31500	1.5	150

Commercial products distributed by UROVISION.  
Not available in all markets!

#### Two-part extension set



Order number	Inner Diameter (mm)	Length (cm)
UD-Z21500	1.5	150

Commercial products distributed by UROVISION.  
Not available in all markets!

#### Infusion tubing



Order number	Length (cm)
UD-Z3000	300

Commercial products distributed by UROVISION.  
Not available in all markets!



**WIRES**

## WIRES

### GUIDE WIRES

#### Standard Guide wires



#### Standard – PTFE coated stainless steel

Order number	Size $\phi$ (inch)	Length (cm)	Tip	Stiffness Tip	Stiffness Shaft
ST-291518	0.018"	150	straight	standard	standard
ST-291525	0.025"	150	straight	standard	standard
ST-291028	0.028"	100	straight	standard	standard
ST-291528	0.028"	150	straight	standard	standard
ST-291532	0.032"	150	straight	standard	standard
ST-291035	0.035"	100	straight	standard	standard
ST-291235	0.035"	125	straight	standard	standard
ST-291535	0.035"	150	straight	standard	stiff
ST-298538	0.038"	85	straight	standard	standard
ST-291038	0.038"	100	straight	standard	standard
ST-291538	0.038"	150	straight	standard	standard
ST-291045	0.045"	100	straight	standard	standard

#### With movable core – PTFE coated stainless steel

Order number	Size $\phi$ (inch)	Length (cm)	Tip	Stiffness Tip	Stiffness Shaft
ST-293501	0.035"	125	straight	standard	standard
ST-293502	0.035"	150	straight	standard	standard
ST-293550	0.035"	150	J-tip	standard	standard

#### Lunderquist – uncoated stainless steel

Order number	Size $\phi$ (inch)	Length (cm)	Tip	Stiffness Tip	Stiffness Shaft
RN-741580	0.035"	80	J-tip (1.0 mm)	standard	standard


#### Lunderquist – PTFE coated stainless steel

Order number	Size $\phi$ (inch)	Length (cm)	Tip	Stiffness Tip	Stiffness Shaft
RN-751580	0.035"	80	J-tip (1.5 mm)	standard	standard
RN-761680	0.035"	80	straight	standard	standard
RN-750380	0.035"	80	J-tip (3.0 mm)	standard	standard
RN-751500	0.035"	100	J-tip (1.5 mm)	standard	standard



## WIRES


### VISIONMULTI



- ◆ Black-and-white coiled
- ◆ Optimal control over advancement
- ◆ Nitinol core increases resistance to kinks
- ◆ PTFE coating enables fast and easy advancement and has optimal gliding properties

Order number	Size $\phi$ (inch)	Length (cm)	Tip	Stiffness Tip	Stiffness Shaft
ST-28150	0.028"	150	straight	flexible	stiff
ST-32150	0.032"	150	straight	flexible	stiff
ST-35150	0.035"	150	straight	flexible	stiff

### VISIONMULTI DUO




- ◆ PTFE coated Nitinol Guide wire with flexibility at both ends
- ◆ Optimally suited for use with flexible ureter endoscopes

Order number	Size $\phi$ (inch)	Length (cm)	Tip	Stiffness Tip	Stiffness Shaft
ST-36150	0.035"	150	straight	flexible	stiff

## WIRES


### VISIOLIX



Attention: Please moisten the wire before use for optimal hydrophilic effect.

- ◆ Hydrophilic coating
- ◆ Super-elastic and stable against kinking
- ◆ Nitinol core

Order number	Size $\varnothing$ (inch)	Length (cm)	Tip	Stiffness Tip	Stiffness Shaft
VLG-67150	0.025"	150	straight	standard	standard
VLG-68150	0.032"	150	straight	standard	standard
VLG-69150	0.035"	150	straight	standard	standard
VLA-71150	0.025"	150	angled	standard	standard
VLA-73150	0.035"	150	angled	standard	standard
VLS-63150	0.035"	150	straight	standard	stiff
VLA-66150-H	0.038"	150	angled	standard	stiff



**ACCESSORIES**

Split cannulas	8.1
Puncture needles	8.2
Dilators	8.3
Adapters	8.4
Diverse	8.7

## ACCESSORIES

# SPLIT CANNULAS

For suprapubic puncture of the bladder

### DESCRIPTION



Order number	For catheter (CH/FR, length)	Length (cm)	Inside diameter (mm)	Outside diameter (mm)
SV-391008	08, 35 cm	8	3.15	3.6
SV-391012	08, 40 cm	12	3.15	3.6
SV-391408	10 +12, 35 cm	8	4.2	4.6
SV-391412	10 +12, 40 cm	12	4.2	4.6
SV-391512	14, 40 cm	12	5.2	5.6

## ACCESSORIES

# NEPHROSTOMY PUNCTURE NEEDLES


### 2-part needle



- ◆ Beveled
- ◆ Sandblasted for optimal ultrasound visibility

Order number	Size (G)	Length (cm)	Maximum Guide wire ø
RN-611820	18	19.5	0.035"


### 2-part needle



- ◆ Trocar tip
- ◆ Sandblasted for optimal ultrasound visibility

Order number	Size (G)	Length (cm)	Maximum Guide wire ø
RN-621820	18	19.5	0.035"
RN-611835	18	35	0.035"


### 2-part Chiba needle



- ◆ Sandblasted for optimal ultrasound visibility

Order number	Size (G)	Length (cm)	Maximum Guide wire ø
RN-632220	22	20	0.018"
RN-632228	22	28	0.018"

### 3-piece needle



- ◆ Beveled
- ◆ Ring markings for optimal ultrasound visibility

Order number	Size (G)	Length (cm)	Maximum Guide wire ø
RN-601820	18	20	0.035"

## ACCESSORIES

# RENOVISION DILATORS



### Dilators / single

Order number	Size (CH/FR)	Length (cm)	Maximum Guide wire ø
RN-770518	05	18	0.035
RN-770618	06	18	0.038"
RN-770718	07	18	0.038"
RN-770818	08	18	0.038"
RN-770918	09	18	0.038"
RN-771018	10	18	0.038"
RN-771218	12	18	0.038"
RN-771418	14	18	0.038"
RN-771618	16	18	0.038"
RN-771818	18	18	0.038"
RN-772018	20	18	0.038"
RN-772218	22	18	0.038"

### Dilator with split shaft / single

Order number	Size (CH/FR)	Length (cm)	Maximum Guide wire ø
RN-780600	06	18	0.038"
RN-780700	07	18	0.038"
RN-780800	08	18	0.038"
RN-780900	09	18	0.038"
RN-781000	10	18	0.038"
RN-781200	12	18	0.038"
RN-781400	14	18	0.038"
RN-781600	16	18	0.038"
RN-781800	18	18	0.038"
RN-782000	20	18	0.038"
RN-782200	22	18	0.038"
RN-782400	24	18	0.038"

### Dilator Set

Order number	Size with split shaft* (CH/FR)	Length (cm)	Maximum Guide wire ø
RN-790816	08/10/12/14/16	18	0.038"
RN-791600	08/10/12/14/16*	18	0.038"

## ACCESSORIES

### ADAPTER

#### UK injection adapter



Order number	Size (G)
AD-270	03
AD-280	04
AD-210	4.7
AD-211	05
AD-212	06
AD-215	07
AD-218	08/09/10

#### Needle adapter with funnel



Order number	Size (G)
AD-310	4.7
AD-311	05
AD-312	06
AD-315	07
AD-318	08/09

#### Tuohy Borst Adapter



Order number	Luer Lock	Suitable for size (CH/FR)
AD-112	Male	06-12
AD-113	Female	06-12

#### Tuohy Borst Adapter



AD-116


Order number	Luer Lock	Suitable for size (CH/FR)
AD-116	With silicone tube (7 cm) and male Luer lock	06-12

#### Universal adapter with step taper and injection apparatus



Order number
AD-103-2

#### Y connector with contrast medium imaging



Order number
RN-550



## ACCESSORIES

### "ONE SIZE FITS ALL" UNIVERSAL ADAPTER

- ◆ Consists of a Luer lock connector and a silicone tube
- ◆ The adapter enables connection of all common urine drainage bag systems with a tapered attachment
- ◆ Especially well suited for connection between Luer lock instruments (nephrostomy and cystostomy)
- ◆ Easy, universal connection possible for patient and medical professional

#### Universal adapter

- ◆ With silicone tube connection, Luer lock male and female
- ◆ With step taper, Luer lock male and female



#### Universal adapter with silicone tube connection, Luer lock male

Order number	Length of connector (cm)	Material property
AD-1005	5	PVC-free
AD-1010	10	PVC-free
AD-1020	20	PVC-free



#### Universal adapter with silicone tube connection, Luer lock female

Order number	Length of connector (cm)	Material property
AD-2005	5	PVC-free



## ACCESSORIES

### ADAPTER

#### Universal adapter with step taper, Luer lock male

Order number	Length of connector (cm)	Material property
AD-100	15	PVC-free
AD-101	-	PVC-free, polyurethane



#### Universal adapter with step taper, Luer lock female

Order number	Length of connector (cm)	Material property
AD-103	-	PVC-free



### Connection adapter

- With step taper and silicone tube connection

#### Connection adapter with step taper and silicone tube connection

Order number	Length of connector (cm)	Material property
AD-220	25	PVC-free



## ACCESSORIES

### DIVERSE

#### Catheter stoppers

##### 070370000

- ◆ Suitable for all balloon catheter funnel attachments

Order number
070370000



Commercial products distributed by UROVISION. Not available in all markets!

#### Urethral olive

Order number
CK-001



Commercial products distributed by UROVISION. Not available in all markets!

#### Block injection

- ◆ Luer injection with a mixture of distilled water and glycerin (10%)

Order number	Capacity (ml)
SP-2100-10	10



#### Wound and bladder injection


- ◆ For wound and bladder lavage, 120 ml

Order number	Capacity (ml)
21120	120



Commercial products distributed by UROVISION. Not available in all markets!

[WWW.UROVISION.DE](http://WWW.UROVISION.DE)



**BIOPSY**


Biopsy needles	9.1
Single-Use Biopsy Gun	9.2
Coaxial needles	9.3

## BIOPSY


# BIOPSY CANNULAS

For the automated removal of biopsied material

DNG-1010 suitable for e.g. PRO-MAG™ 2.2



DNG-1020 suitable for e.g. PRO-MAG™ Ultra



- ◆ First-class beveling process for high-quality histological findings
- ◆ Centimeter gradation for exact determination of incision depth
- ◆ Echoic tip for clear visibility in ultrasound
- ◆ Size easily identified based on color marking

Order number	Size (G)	Length (cm)	Suitable Coaxial needle
DNG-1010-14-0100	14	10	COX-1010-14-0046
DNG-1010-14-0120	14	12	
DNG-1010-14-0160	14	16	COX-1010-14-0106
DNG-1010-14-0200	14	20	
DNG-1010-16-0100	16	10	COX-1010-16-0046
DNG-1010-16-0160	16	16	COX-1010-16-0106
DNG-1010-16-0200	16	20	
DNG-1010-16-0250	16	25	
DNG-1010-18-0100	18	10	COX-1010-18-0049
DNG-1010-18-0160	18	16	COX-1010-18-0106
DNG-1010-18-0200	18	20	COX-1010-18-0146
DNG-1010-18-0250	18	25	
DNG-1010-18-0300	18	30	
DNG-1010-20-0100	20	10	COX-1010-20-0049
DNG-1010-20-0160	20	16	COX-1010-20-0109
DNG-1010-20-0200	20	20	

Order number	Size (G)	Length (cm)	Suitable Coaxial needle
DNG-1020-14-0100	14	10	COX-1020-14-0046
DNG-1020-14-0120	14	12	
DNG-1020-14-0160	14	16	COX-1020-14-0106
DNG-1020-14-0200	14	20	
DNG-1020-16-0100	16	10	COX-1020-16-0046
DNG-1020-16-0160	16	16	COX-1020-16-0106
DNG-1020-16-0200	16	20	
DNG-1020-16-0250	16	25	
DNG-1020-18-0100	18	10	COX-1020-18-0046
DNG-1020-18-0160	18	16	COX-1020-18-0106
DNG-1020-18-0200	18	20	COX-1020-18-0146
DNG-1020-18-0250	18	25	
DNG-1020-18-0300	18	30	
DNG-1020-20-0100	20	10	COX-1020-20-0049
DNG-1020-20-0160	20	16	COX-1020-20-0109
DNG-1020-20-0200	20	20	


Commercial products distributed by UROVISION. Not available in all markets!

## BIOPSY

# BIOPSIE CANNULAS

For the automated removal of biopsied material

DNG-1040 suitable for e.g. Bard® Magnum®



- ◆ First-class beveling process for high-quality histological findings
- ◆ Centimeter gradation for exact determination of incision depth
- ◆ Echoic tip for clear visibility in ultrasound
- ◆ Size easily identified based on color marking

Order number	Size (G)	Length (cm)	Suitable Coaxial needle
DNG-1040-14-0100	14	10	COX-1040-14-0043
DNG-1040-14-0130	14	13	
DNG-1040-14-0160	14	16	COX-1040-14-0103
DNG-1040-14-0200	14	20	
DNG-1040-16-0100	16	10	COX-1040-16-0043
DNG-1040-16-0130	16	13	
DNG-1040-16-0160	16	16	COX-1040-16-0103
DNG-1040-16-0200	16	20	
DNG-1040-18-0100	18	10	COX-1040-18-0043
DNG-1040-18-0130	18	13	
DNG-1040-18-0160	18	16	COX-1040-18-0103
DNG-1040-18-0200	18	20	COX-1040-18-0143
DNG-1040-18-0250	18	25	
DNG-1040-18-0300	18	30	
DNG-1040-20-0100	20	10	COX-1040-20-0046
DNG-1040-20-0130	20	13	
DNG-1040-20-0160	20	16	COX-1040-20-0106
DNG-1040-20-0200	20	20	

Commercial products distributed by UROVISION. Not available in all markets!




## BIOPSY

# DISPOSABLE BIOPSY GUN

### For automatic biopsy

#### ABG-2020



- ◆ Optimized cut to obtain high quality histological samples
- ◆ Centimeter markings for exact positioning
- ◆ Echogenic tip for clear visibility using ultrasound
- ◆ One-hand-use with choice of driver either on side or back
- ◆ Can be used semi- or full automatic
- ◆ Side-notch 19 mm

Order number	Size (G)	Length (cm)	suitable for coaxial Needle
ABG-2020-16-0100	16	10	COX-2000-16-0051
ABG-2020-16-0150	16	15	COX-2000-16-0101
ABG-2020-16-0200	16	20	COX-2000-16-0151
ABG-2020-18-0100	18	10	COX-2000-18-0051
ABG-2020-18-0150	18	15	COX-2000-18-0101
ABG-2020-18-0200	18	20	COX-2000-18-0151
ABG-2020-18-0250	18	25	COX-2000-18-0201


Commercial products distributed by UROVISION. Not available in all markets!

## BIOPSY

# DISPOSABLE BIOPSY GUN

For automatic biopsy

ABG-2040



- ◆ Optimized cut to obtain high quality histological samples
- ◆ Centimeter markings for exact positioning
- ◆ Echogenic tip for clear visibility using ultrasound
- ◆ One-hand-use with choice of driver either on side or back
- ◆ full automatic
- ◆ Side-notch 20 mm
- ◆ Strong spring tension
- ◆ Greater patient comfort thanks to special low-noise attenuation

Order number	Size (G)	Length (cm)	suitable for coaxial Needle
ABG-2040-16-0100	16	10	COX-2000-16-0051
ABG-2040-16-0150	16	15	COX-2000-16-0101
ABG-2040-16-0200	16	20	COX-2000-16-0151
ABG-2040-18-0100	18	10	COX-2000-18-0051
ABG-2040-18-0150	18	15	COX-2000-18-0101
ABG-2040-18-0200	18	20	COX-2000-18-0151
ABG-2040-18-0250	18	25	COX-2000-18-0201

Commercial products distributed by UROVISION. Not available in all markets!

## BIOPSY

# COAXIAL NEEDLE

- ◆ Allows multiple sampling through only one percutaneous puncture
- ◆ Centimeter graduation and infinitely variable spring clips for exact determination of penetration depth
- ◆ Echogenic tip for clear ultrasound visibility
- ◆ Luer Lock adapter for anesthetic application

## COX



Order number	Size (CH/FR)	Length (cm)
COX-1010-14-0046	13	4,6
COX-1010-14-0106	13	10,6
COX-1010-16-0046	15	4,6
COX-1010-16-0106	15	10,6
COX-1010-18-0049	17	4,6
COX-1010-18-0106	17	10,6
COX-1010-18-0146	17	14,6
COX-1010-20-0049	19	4,9
COX-1010-20-0109	19	10,9
COX-1020-14-0046	13	4,6
COX-1020-14-0106	13	10,6
COX-1020-16-0046	15	4,6
COX-1020-16-0106	15	10,6
COX-1020-18-0046	17	4,6
COX-1020-18-0106	17	10,6
COX-1020-18-0146	17	14,6
COX-1020-20-0049	19	4,9
COX-1020-20-0109	19	10,9
COX-1040-14-0043	13	4,3
COX-1040-14-0103	13	10,3
COX-1040-16-0043	15	4,3
COX-1040-16-0103	15	10,3
COX-1040-18-0043	17	4,3
COX-1040-18-0103	17	10,3
COX-1040-18-0143	17	14,3
COX-1040-20-0046	19	4,6
COX-1040-20-0106	19	10,6

Order number	Size (CH/FR)	Length (cm)
COX-2000-14-0051	13	5,1
COX-2000-14-0101	13	10,1
COX-2000-14-0151	13	15,1
COX-2000-16-0051	15	5,1
COX-2000-16-0101	15	10,1
COX-2000-16-0151	15	15,1
COX-2000-18-0051	17	5,1
COX-2000-18-0101	17	10,1
COX-2000-18-0151	17	15,1
COX-2000-18-0201	17	20,1
COX-2000-20-0051	19	5,1
COX-2000-20-0101	19	10,1
COX-2000-20-0151	19	15,1

Commercial products distributed by UROVISION.  
 Not available in all markets!

[WWW.UROVISION.DE](http://WWW.UROVISION.DE)




**SUPRAPUBIC  
BLADDER DRAINAGE**

SUPRAVISION PUR Set	10.1
SUPRAVISION TUR Puncture Set	10.2
SUPRAVISION TUR Pigtail Catheter	10.3
SUPRAVISION Integral Standard Suprapubic catheter with integrated balloon	10.4
SUPRAVISION Balloon Catheter & Sets	10.5

## SUPRAPUBIC BLADDER DRAINAGE

# SUPRAVISION PUR Puncture and Exchange Set with pigtail and funnel

### SET COMPONENTS



### Pigtail Catheter

- ◆ Polyurethane
- ◆ Flexible pigtail prevents irritation
- ◆ Funnel for attaching a urine drainage bag
- ◆ Markings for exact positioning




Puncture Set				
Order number	Size (CH/FR)	Length (cm)	Split cannula diameter size (mm)	Split cannula length (mm)
SV-378565	08	65	3.6	120
SV-371065	10	65	4.6	120
SV-371265	12	65	5.6	120

Exchange Set						
Order number	Size (CH/FR)	Length (cm)	Wire ø (inch)	Wire length (cm)	Wire tip	Wire stiffness
SV-371066	10	65	0.038"	100	straight	standard
SV-371266	12	65	0.038"	100	straight	standard
SV-371466	14	65	0.038"	100	straight	standard

## SUPRAPUBIC BLADDER DRAINAGE

### SUPRAVISION TUR Puncture Set with pigtail and funnel

#### SET COMPONENTS



#### Pigtail Catheter

- ◆ Polyurethane
- ◆ Open or closed tip
- ◆ Optimal for draining and rinsing during transurethral resection
- ◆ Flexible pigtail prevents irritation
- ◆ Funnel for attaching a urine drainage bag
- ◆ Markings for exact positioning



Puncture Set					
Order number Open tip	Order number Closed tip	Size (CH/FR)	Length (cm)	Split cannula diameter size (mm)	Split cannula length (mm)
SV-391540	SV-391640	14	65	5.6	120



## SUPRAPUBIC BLADDER DRAINAGE

# SUPRAVISION TUR PIGTAIL CATHETER


With pigtail and funnel

### DESCRIPTION



### TUR Pigtail Catheter

- ◆ Polyurethane
- ◆ Separate funnel for attachment
- ◆ Open tip
- ◆ Optimal for draining and rinsing during TUR
- ◆ Flexible pigtail prevents irritation
- ◆ Excellent urine routing due to extra-large drainage openings in the pigtail
- ◆ Detachable funnel for connecting a urine drainage bag
- ◆ Markings for exact positioning
- ◆ Suitable for reusable trocars




Order number	Size (CH/FR)	Length (cm)
SV-391865-ML	14	65

## SUPRAPUBIC BLADDER DRAINAGE

# SUPRAVISION INTEGRAL STANDARD

Suprapubic catheter with integrated balloon, 100% silicone

### SET COMPONENTS



### Suprapubic Catheter

- ◆ Completely silicone catheter with integrated balloon, enabling easier catheter placement
- ◆ 100% silicone, transparent with an open tip
- ◆ Excellent flow-through rate
- ◆ Very low propensity for encrustation
- ◆ Low-risk, easy handling
- ◆ Rounded tip for easier catheter placement
- ◆ Reinforced funnel attachment to prevent funnel from ripping and to enable secure connection of the urine drainage bag
- ◆ Transparent with integrated x-ray contrast strips
- ◆ Very good lubricating properties
- ◆ Stable inner and outer coating



### VISIOSLIDE Special Coating

- ◆ Improved lubricating properties, wear-resistant inside and outside
- ◆ Lower propensity for encrustation, resulting in increased life cycle

## SUPRAPUBIC BLADDER DRAINAGE

# SUPRAVISION INTEGRAL STANDARD

Suprapubic catheter with integrated balloon, 100% silicone

Single catheter					
Order number Normal variant	Order number Normal variant with VISIOSLIDE	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Suitable for wire ø (inch)
SV-310835	SV-310835-SC	8	35	1.5	0.028"
SV-311040	SV-311040-SC	10	40	3	0.032"
SV-311240	SV-311240-SC	12	40	5	0.038"
SV-311440	SV-311440-SC	14	40	5	0.038"
SV-311640	SV-311640-SC	16	40	10	0.038"
SV-311840	SV-311840-SC	18	40	10	0.038"
SV-312040	SV-312040-SC	20	40	10	0.038"
SV-312240	SV-312240-SC	22	40	10	0.038"
SV-312440	SV-312440-SC	24	40	10	0.038"

Single catheter					
Order number Short variant	Order number Short variant with VISIOSLIDE	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Suitable for wire ø (inch)
SV-321025	SV-321025-SC	10	25	3	0.032"
SV-321225	SV-321225-SC	12	25	5	0.038"
SV-321425	SV-321425-SC	14	25	5	0.038"
SV-321625	SV-321625-SC	16	25	10	0.038"
SV-321825	SV-321825-SC	18	25	10	0.038"
SV-322025	SV-322025-SC	20	25	10	0.038"
SV-322225	SV-322225-SC	22	25	10	0.038"

Bladder Puncture Set				
Order number With scalpel	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Wire ø (inch)
SV-300835	8	35	1.5	0.028"
SV-301025	10	25	3	0.032"
SV-301040	10	40	3	0.032"
SV-301225	12	25	5	0.038"
SV-301240	12	40	5	0.038"
SV-301440	14	40	5	0.038"

## SUPRAPUBIC BLADDER DRAINAGE

### SUPRAVISION INTEGRAL STANDARD

Suprapubic catheter with integrated balloon, 100% silicone

Initial Bladder Puncture Set with VISIOSLIDE Special Coating				
Order number With scalpel With VISIOSLIDE	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Wire $\phi$ (inch)
SV-300835-SC	8	35	1.5	0.028"
SV-301040-SC	10	40	3	0.032"
SV-301240-SC	12	40	5	0.038"
SV-301440-SC	14	40	5	0.038"

Exchange Set					
Order number Normal variant	Order number Normal variant with VISIOSLIDE	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Wire $\phi$ (inch)
SV-330835	SV-330835-SC	8	35	1.5	0.028"
SV-331040	SV-331040-SC	10	40	3	0.032"
SV-331240	SV-331240-SC	12	40	5	0.038"
SV-331440	SV-331440-SC	14	40	5	0.038"
SV-331640	SV-331640-SC	16	40	10	0.038"
SV-331840	SV-331840-SC	18	40	10	0.038"
SV-332040	SV-332040-SC	20	40	10	0.038"


Dilation exchange set				
Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Wire $\phi$ (inch)
SV-341040	10	40	3	0.032"
SV-341240	12	40	5	0.038"
SV-341440	14	40	5	0.038"
SV-341640	16	40	10	0.038"
SV-341840	18	40	10	0.038"

## SUPRAPUBIC BLADDER DRAINAGE

# SUPRAVISION SOFT

Silicone-coated balloon catheter with latex core

### SET COMPONENTS



### Balloon Catheter

- ◆ The advantages of latex and silicone are combined in one catheter by using a silicone layer on a latex core
- ◆ The funnel attachment does not fall off and is optimally secured on the bag system connector
- ◆ Soft balloon catheter with a high level of wearing comfort for the patient, so it is ideally suited for long-term drainage
- ◆ The catheter has very high balloon stability, largely eliminating the possibility of a rupture or subsequent blockage
- ◆ The highly rounded tip enables easy placement of catheter
- ◆ Open tip




Order number Single catheter	Order number Exchange set	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Wire ø (inch)
SV-361240	SV-351240	12	40	10	0.038"
SV-361440	SV-351440	14	40	10	0.038"
SV-361640	SV-351640	16	40	10	0.038"
SV-361840	SV-351840	18	40	10	0.038"
SV-362040	SV-352040	20	40	10	0.038"
SV-362240	-	22	40	10	0.038"
SV-362440	-	24	40	10	0.038"
SV-362640	-	26	40	10	0.038"

## SUPRAPUBIC BLADDER DRAINAGE

### SUPRAVISION INTEGRAL PIGTAIL

Balloon catheter with pigtail tip, 100% silicone

#### SET COMPONENTS



#### Balloon Catheter

- ◆ 6 drainage openings in the pigtail, providing very high drainage performance
- ◆ Silicone for optimal tissue compatibility
- ◆ Less irritation when placing and exchanging the catheters due to integrated balloon
- ◆ Very low propensity for encrustation
- ◆ Reinforced funnel attachment to prevent funnel from ripping and to enable secure connection of the urine drainage bag
- ◆ With integrated balloon
- ◆ 100% silicone



## SUPRAPUBIC BLADDER DRAINAGE

# SUPRAVISION INTEGRAL PIGTAIL

### Balloon catheter with pigtail tip, 100% silicone

Single catheter				
Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Suitable for wire ø (inch)
SV-401240	12	40	5	0.038"
SV-401440	14	40	5	0.038"
SV-401640	16	40	10	0.038"
SV-401840	18	40	10	0.038"

Puncture Set				
Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Suitable for wire ø (inch)
SV-411240	12	40	5	0.038"
SV-411440	14	40	5	0.038"


Exchange Set				
Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)	Wire ø (inch)
SV-421240	12	40	5	0.038"
SV-421440	14	40	5	0.038"
SV-421640	16	40	10	0.038"

## SUPRAPUBIC BLADDER DRAINAGE

### SUPRAVISION INTEGRAL SPECIAL TIP

Suprapubic dilatation balloon catheter with flat integrated balloon, 100% silicone

#### DESCRIPTION



#### Dilatation balloon catheter

- ◆ Completely 100% silicone catheter with flat integrated balloon, enabling easier catheter placement
- ◆ Shaped and reinforced open tip
- ◆ Dilatation is facilitated by the metal-reinforced tip and exact adaptation to the Guide wire
- ◆ Very low incidence of encrustation, low-risk, and easy application
- ◆ Optional stylet; possible to probe the puncture site again if catheter is dislocated. This can help to avoid repeated punctures

Order number Without stylet	Order number With stylet	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
SV-441240	SV-511240	12	40	5
SV-441440	SV-511440	14	40	5
SV-441640	SV-511640	16	40	10
SV-441840	SV-511840	18	40	10
SV-442040	SV-512040	20	40	10

Suitable for 0.035" Guide wire Wire order number: ST-291035





**TRANSURETHRAL  
CATHETERS**


Transurethral Catheters	11.1
Extension Catheters	11.7
VISIOKON Catheter	11.8

## TRANSURETHRAL CATHETERS

# VISI-ESOFT BALLOON CATHETER

Extra soft, 100% silicone

### DESCRIPTION



### Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Very soft balloon catheter for low irritation and a very high comfort level for the patient
- ◆ Very low propensity for encrustation
- ◆ Very good long-term draining performance due to large inner diameter


Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KG-300810	8	30	5
KG-301010	10	30	5
KG-301210	12	40	10
KG-301410	14	40	10
KG-301610	16	40	10
KG-301810	18	40	10
KG-302010	20	40	10
KG-302210	22	40	10
KG-302410	24	40	10

## TRANSURETHRAL CATHETER

# VISI-INTEGRAL BALLOON CATHETER

With integrated balloon, 100% silicone

## DESCRIPTION



## Balloon Catheter

- ◆ Gently dilating Nelaton tip
- ◆ Silicone for optimal tissue compatibility
- ◆ Less irritation in the urethra when placing and exchanging the catheters due to integrated balloon
- ◆ Very low propensity for encrustation
- ◆ Reinforced funnel attachment to prevent funnel from ripping and to enable secure connection of the urine drainage bag


Order number Nelaton tip, closed	Order number Nelaton tip with end hole	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-101040	KI-101040-SL	10	30	3
KI-101240	KI-101240-SL	12	40	5
KI-101440	KI-101440-SL	14	40	10
KI-101640	KI-101640-SL	16	40	10
KI-101840	KI-101840-SL	18	40	10
KI-102040	KI-102040-SL	20	40	10
KI-102240	KI-102240-SL	22	40	10
KI-102440	KI-102440-SL	24	40	10
KI-102640	KI-102640-SL	26	40	10

## TRANSURETHRAL CATHETER

# VISISOFT BALLOON CATHETER

With latex core and solid silicone coating inside and outside

### DESCRIPTION



### Balloon Catheter

- ◆ Nelaton tip
- ◆ Silicone coating for optimal tissue compatibility
- ◆ High balloon stability, which largely eliminates the possibility of balloon ruptures and subsequent blockage
- ◆ Soft balloon catheter offers wearing comfort for patients and is well suited for long-term drainage
- ◆ Non-rip funnel attachment offers ideal anchoring for bag system step connector
- ◆ Bladder stone resistance of the balloon


Order number Nelaton tip, closed	Order number Nelaton tip with end hole	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-200810	-	8	30	5
KI-201010	-	10	30	5
KI-201210	KI-201210-SL	12	40	10
KI-201410	KI-201410-SL	14	40	10
KI-201610	KI-201610-SL	16	40	10
KI-201810	KI-201810-SL	18	40	10
KI-202010	KI-202010-SL	20	40	10
KI-202210	KI-202210-SL	22	40	10
KI-202410	KI-202410-SL	24	40	10
KI-202610	KI-202610-SL	26	40	10

## TRANSURETHRAL CATHETER

# VISI-INTEGRAL SOFT BALLOON CATHETER

100% silicone

## DESCRIPTION



## Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Secure placement of the catheter in the bladder due to special balloon properties
- ◆ Easier placement and exchange of catheters due to integrated balloon
- ◆ Very low propensity for encrustation
- ◆ Balloon catheter with rotation-secure atraumatic Tiemann tip
- ◆ Optimally suited for long and short-term drainage


Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KI-401010	10	31	5
KI-401210	12	40	5
KI-401410	14	40	5
KI-401610	16	40	10
KI-401810	18	40	10
KI-402010	20	40	10
KI-402210	22	40	10
KI-402410	24	40	10

## TRANSURETHRAL CATHETER

# VISIOCATH BALLOON CATHETER

100% silicone

## DESCRIPTION



### Balloon Catheter

- ◆ With VISIOSLIDE special coating or uncoated (standard)
- ◆ Nelaton tip
- ◆ For women, men, and children
- ◆ Silicone for optimal tissue compatibility
- ◆ High balloon stability with low diffusion
- ◆ Reinforced funnel attachment to prevent funnel from ripping and to enable secure connection of the urine drainage bag
- ◆ Transparent with integrated x-ray contrast strips
- ◆ Also available as a children's catheter with stiffening mandrin
- ◆ Very good lubricating properties
- ◆ Stable inner and outer coating

### VISIOSLIDE Special Coating

- ◆ Improved lubricating properties, wear-resistant inside and outside
- ◆ Lower propensity for encrustation, resulting in increased life cycle

## TRANSURETHRAL CATHETER

# VISIOCATH BALLOON CATHETER

100% silicone

Balloon Catheter for Men – Standard				
Order number Nelaton tip, closed	Order number Nelaton tip with end hole	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-441210	KW-441210-SL	12	40	10
KW-441410	KW-441410-SL	14	40	10
KW-441610	KW-441610-SL	16	40	10
KW-441810	KW-441810-SL	18	40	10
KW-442010	KW-442010-SL	20	40	10
KW-442210	KW-442210-SL	22	40	10
KW-442410		24	40	10
KW-442610		26	40	10
KW-441430	KW-441430-SL	14	40	30
KW-441630	KW-441630-SL	16	40	30
KW-441830	KW-441830-SL	18	40	30
KW-442030	KW-442030-SL	20	40	30
KW-442230	KW-442230-SL	22	40	30
KW-442430	KW-442430-SL	24	40	30

Balloon Catheter for Men with Visioslide Special Coating			
Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-441210-SC	12	40	10
KW-441410-SC	14	40	10
KW-441610-SC	16	40	10
KW-441810-SC	18	40	10
KW-442010-SC	20	40	10
KW-442210-SC	22	40	10
KW-442410-SC	24	40	10



## TRANSURETHRAL CATHETER

# VISIOCATH BALLOON CATHETER

100% silicone


Balloon Catheter for Women			
Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KW-451210	12	22	10
KW-451410	14	22	10
KW-451610	16	22	10
KW-451810	18	22	10
KW-452010	20	22	10
KW-452210	22	22	10
KW-452410	24	22	10

## TRANSURETHRAL CATHETER

# WHISTLE TIP CATHETER

100% silicone

### DESCRIPTION



### Whistle tip catheter


- ◆ With atraumatic, open and short whistle tip
- ◆ Especially for patients receiving enteral feeding
- ◆ Silicone for optimal tissue compatibility
- ◆ Reliable drainage performance due to large inside diameter for problematic urine findings

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
KF-501610	16	40	10
KF-501810	18	40	10
KF-502010	20	40	10

## TRANSURETHRAL CATHETER

# UROVISION EXTENSION CATHETER

### DESCRIPTION



### Extension catheter

- ◆ Nongraduated expansion of urethra, so there is only a single lesion in the mucous membrane
- ◆ Patient-friendly longitudinal passage (pain-free)
- ◆ Low recurrence rate of urethral strictures

Straight, open tip		
Order number	Size (CH/FR)	Length (cm)
EK-811022	10-22	40

Tiemann tip		
Order number	Size (CH/FR)	Length (cm)
EK-801022	10-22	40

## TRANSURETHRAL CATHETER

### VISIOKON CATHETER

100% silicone, for x-ray imaging of female urethra

#### DESCRIPTION



#### Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ X-ray contrast strips provide optimal visibility along the entire length of the catheter
- ◆ Pinch clamp for securing the cone
- ◆ Catheter and cone made of 100% silicone with additional Visioslide coating for improved lubricating properties
- ◆ Optimum contrast medium visualization due to large eyes at catheter end

Order number	Size (CH/FR)	Length (cm)	Balloon capacity (ml)
PI-230	14	40	10



APPROXIMATE ML  
2000  
1800  
1600  
1400

**URELIA®**  
by Urovision

F3

100  
75  
50  
25

**URINE DRAINAGE  
SYSTEMS**

Urine drainage bags	12.1
VELCRO Straps	12.6
Drainage bag Holders & One-Leg Pants	12.7

## URINE DRAINAGE SYSTEMS

# COMFORT BEDSIDE BAG

## Closed Urine Drainage System and Bedside Bag

### DESCRIPTION



### Comfort Bedside Bag

- ◆ Gradation and odor-proof bag foil
- ◆ Drip chamber with foil anti-reflux valve provides optimal protection from rising germs
- ◆ One-hand tap with open holder
- ◆ Stable, removable nylon hook and cord for fastening to the bed
- ◆ Sliding clamp
- ◆ Bed sheet clip
- ◆ Soft, kink-resistant tube
- ◆ Urine sampling port
- ◆ Sterilized

Comfort Bedside Bag		
Order number	Bag volume (ml)	Tube attachment (cm)
UVD-47587	2000	120

## URINE DRAINAGE SYSTEMS

# URELIA – CLOSED URINE DRAINAGE SYSTEM

Urine drainage system for long-term care

## DESCRIPTION



### UVD-23

- ◆ 2,000 ml volume
- ◆ Hose length 120 cm, kink resistant with sliding clamp
- ◆ Needle-free urine sampling site
- ◆ Drip chamber with integrated non-return valve
- ◆ Bacteria-proof discharge valve
- ◆ Simple one-hand drainage
- ◆ Double hook attachment
- ◆ Bed sheet clamp
- ◆ Sterile

### Closed urine drainage system

Order number	Bag volume (ml)	Hose length (cm)
UVD-23	2000	120



## URINE DRAINAGE SYSTEMS

# URELIA – URINE DRAINAGE SYSTEM

Urine drainage system for short-term use

## DESCRIPTION



### UVD-43

- ◆ 2,000 ml volume
- ◆ Hose length 110 cm, kink resistant with sliding clamp
- ◆ Needle-free urine sampling site
- ◆ Return barrier
- ◆ Bacteria-proof discharge valve
- ◆ Simple one-hand drainage
- ◆ Conical step connector for optimum hold
- ◆ Sterile



### plastic hanger

Order number

M02175NS

### Urine drainage system for short-term use

Order number	Bag volume (ml)	Hose length (cm)
UVD-43	2000	110

## URELIA – DAY / NIGHT BAGS

### UVD-1030/1060/1090

#### DESCRIPTION



#### The secure all-round care for day and night

- ◆ Sturdy day/night bags with 1,000 ml volume
- ◆ Special, skin-friendly, soft back of the bag
- ◆ Coated hose to avoid irritation
- ◆ Different hose lengths available
- ◆ Needle-free urine sampling site
- ◆ Drip chamber with integrated non-return valve
- ◆ Bacteria-proof discharge valve
- ◆ Dog house for easy storage of the drain tap
- ◆ Simple one-hand drainage
- ◆ Sterile


Order number	Bag volume (ml)	Hose length (cm)
UVD-1030	1000	30
UVD-1060	1000	60
UVD-1090	1000	90

## URINE DRAINAGE SYSTEMS

# URELIA LEG BAG WITH DOG HOUSE

## UVD-50XX / UVD-75XX

### DESCRIPTION



### Wide range of high-quality, excellent leg bags with dog house

- ◆ Stable leg bags with 500 ml or 750 ml volume
- ◆ Hose length 10 cm or 50 cm (can be shortened with exchangeable universal adapter)
- ◆ Coated hose to avoid irritation
- ◆ Special, skin-friendly, soft back of the bag
- ◆ Dog house for easy storage of the drain tap
- ◆ Washable special fastening straps for secure attachment on the leg (available separately)
- ◆ Can be connected to bed bag
- ◆ Simple one-hand drainage
- ◆ Return valve
- ◆ Sterile

Order number	Bag volume (ml)	Hose length (cm)
UVD-5010	500	10
UVD-5050	500	50 (toshorten)
UVD-7510	750	10
UVD-7550	750	50 (toshorten)

## URELIA – 3-CHAMBER LEG BAG

for high wearing comfort with inconspicuous fit

### DESCRIPTION



#### 3-chamber leg bag

- ◆ Stable urine bags with 500 ml or 750 ml volume
- ◆ Adapts inconspicuously to the shape of the leg thanks to the three chambers
- ◆ Special, skin-friendly, soft back of the bag
- ◆ Coated hose to avoid irritation
- ◆ Hose length 50 cm, can be shortened with adapter
- ◆ Simple one-hand drainage
- ◆ Return valve
- ◆ Can be connected to bed bag
- ◆ Washable special fastening straps for secure attachment to the leg (available separately)
- ◆ Sterile

Order number	Bag volume (ml)	Hose length (cm)
UVD-3K50	500	50 (kürzbar)
UVD-3K75	750	50 (kürzbar)

## FASTENING STRAPS FOR LEG BAGS (AVAILABLE SEPARATELY)

For attachment of leg bags

Order number
UVD-552500



- ◆ Set of two straps
- ◆ 5.5 cm wide strap and 2.5 cm wide strap
- ◆ With silicone coating inside for optimum grip
- ◆ Washable at 30° Celsius

Order number
UVD-591005



- ◆ Elastic universal strap with non-slip silicone strips and
- ◆ Velcro fastener (thigh attachment)
- ◆ Basic length 50 cm, width 2.5 cm
- ◆ (VE: 1 pair)

## URINE DRAINAGE SYSTEMS

# BEDSIDE BAG HOLDERS & ONE-LEG PANTS

### Bedside Bag Holders

- ◆ Bedside bag holders for all bedside bag sizes and shapes
- ◆ Can be closed with velcro strap
- ◆ Discreet case for the urine drainage bag




Order number
BETTBLTA FEMALE
BETTBLTA MALE
BETTBLTA NEUTRAL

### One-Leg Pant for Leg Bag

- ◆ Leg bag holder for commercially available leg bags
- ◆ Can be washed at 95°C
- ◆ Skin-friendly and stretchable waistband
- ◆ High level of wearing comfort, with no constriction of leg
- ◆ Opening for tap for easier emptying of the bag
- ◆ Secure fastening of the urine bag
- ◆ Over-the-counter aid

Order number	Size	Waist (cm)
EB-15611	S (severity level)	50 - 64
EB-15612	M	65 - 85
EB-15613	L	86 - 100
EB-15614	XL	Over 101



Commercial products distributed by UROVISION. Not available in all markets!

[WWW.UROVISION.DE](http://WWW.UROVISION.DE)



**UROVISION<sup>®</sup>**  
A MEDI-GLOBE COMPANY

Urine Drainage  
Systems

Transurethral  
Catheters

Suprapubic Bladder  
Drainage

Biopsy

Accessories

Wires

Urodynamics

Surgical Catheters

Ureteral Catheters

Stone Baskets

Nephrostomy

Ureteral Stents



# We are there for you around the globe

UROVISION belongs to the Medi-Globe Group, a company with high commitment, great efficiency and a solid financial structure. With 150 employees at our headquarters in Achenmühle near Rosenheim, our own branches in the Czech Republic and Austria and worldwide trade partnerships, we ensure that our modern medical products for urology are available to doctors, healthcare professionals and, of course, patients around the globe.

We are learning every day. The changing legal requirements in the countries, the evolving requirements, the medical or medical technological advancements: Thanks to our worldwide network, we are always up to date – and able to react quickly.

Cooperation based on partnership is one of our great strengths. We are open-minded for all players in our field of activity.

We are URO.  
Simply give us a call!

## **UROVISION GmbH**

Medi-Globe-Straße 1-5  
D-83101 Rohrdorf OT Achenmühle  
Phone: +49 (0) 8032 973-400  
Fax: +49 (0) 8032 973-411  
E-Mail: [info@urovision.com](mailto:info@urovision.com)  
[www.urovision-urotech.de](http://www.urovision-urotech.de)