

# PRIMUS®

Advanced Single-Use Endoscopic Imaging  
with PRIMUS® Uretero and PRIMUS® Cysto




- ✓ Consistently sterile
- ✓ Consistently safe
- ✓ Consistently suitable

*Your Partner in the  
Endourology*

# PRIMUS® Single-use flexible Ureteroscope

Single-use ureteroscope with 7.5 CH/FR outer diameter and included suction function

- ✓ Innovative
- ✓ Precise
- ✓ Versatile




- ✓ Bending lock to lock deflection of distal tip
- ✓ 1:1 torque
- ✓ Rotatable sleeve for hard to reach areas and increased precision
- ✓ Ergonomic handle
- ✓ Integrated suction function
- ✓ Fixated Tuohy Borst adapter

Article-no.	Specifications				
	Outer diameter tube	Working length of shaft	Size working channel	Deflection (up/down)	Field of view
URS-7.5S	7.5 CH/FR	700 mm	3.6 CH/FR	Standard deflection 285°/285°	120°
URS-7.5S-R				Reverse deflection 285°/285°	

# PRIMUS® Single-use flexible Cystoscope

Single-use high definition Cystoscope with 15.3 CH/FR outer diameter

- ✓ Efficient
- ✓ Reliable
- ✓ Easy handling



- ✓ 210 degree deflection in both directions
- ✓ Included Tuohy Borst adapter
- ✓ HD imaging quality
- ✓ Ergonomic handle

Article-no.	Specifications				
	Outer diameter tube	Working length of shaft	Size working channel	Deflection (up/down)	Field of view
CYSHD-16	15.3 CH/FR	400 mm	6.6 CH/FR	Standard deflection 210°/210°	110°
CYSHD-16-R				Reverse deflection 210°/210°	

# PRIMUS® Monitor

Digital video HD-Monitor for all PRIMUS® endoscopic devices



- ✓ Integrated battery function
- ✓ Compatible with all PRIMUS® scopes
- ✓ DICOM connectivity
- ✓ Allows usage of two scopes simultaneously
- ✓ Use as stand-alone monitor or in combination with your video tower

Article-no.	Specifications				
	Max. resolution	Display	Internal storage capacity	Digital video output	Integrated battery
MEDVM-D1	1920x1080	15.6" TFT LCD	64 GB	HDMI and SDI	4 hours

#### Manufacturer

Hunan Vathin Medical Instrument Co.  
1/F, Building 12,  
Innovation and Entrepreneurship  
Service Center,  
No 9 Chuanqi West Road, Jiuhua  
Economic Development Zone,  
Xiangtan, Hunan, China

#### Exclusive Distributor

Urotech GmbH  
Medi-Globe-Str. 1-5  
DE-83101 Rohrdorf OT Achenmühle  
Germany  
Telephone: +49 (0) 8032 - 973-210  
Fax: +49 (0) 8032 - 973-211  
E-Mail: [info@urotech.com](mailto:info@urotech.com)  
[www.urotech.com](http://www.urotech.com)