

UROCKIT® Steering Device

Single-use steering device for the controlled placement of ureteral stents

- ✓ Compatible with Urotech ureteral stents
- ✓ For guidewires up to max. 0.035"
- ✓ Radiopaque marker
- ✓ 42 cm working length
- ✓ Sterile single use

1 Advance

2 Rotate

3 Release

UROCKIT® Steering Device

The UROCKIT® Steering Device supports the controlled advancement, positioning and final release of compatible ureteral stents during implantation.


- ✓ **Controlled placement** of ureteral stents
- ✓ **Supports precise positioning** in the renal pelvis and urinary bladder
- ✓ **Stabilises the ureteral stent** during guidewire removal
- ✓ **Supports position verification** using imaging
- ✓ **Allows correction of the ureteral stent position** during implantation
- ✓ **Supports reliable urinary drainage** through correct ureteral stent positioning

Technical specifications

FEATURE	SPECIFICATION
Product type	Sterile single-use steering device
Working length	42 cm
Guidewire	max. 0.035"
Radiopacity	M variants with marker
Shelf life	3 years
Packaging unit	10 pieces
eIFU	available
MDR class	Is


Compatibility with Urotech products

SD-0750-M: compatible with 6, 7 und 8 CH/FR ureteral stents
SD-0550-M: compatible with 5 CH/FR ureteral stents



UROCKIT® Steering Device

For a controlled connection between the ureteral stent and the steering device. Suitable for all Urotech ureteral stents from 5-8 CH/FR.



Dual-tube system

Supports advancement, rotation and release

One-handed operation

For easy release of the ureteral stent

Handle design

For ergonomic handling

Defined connection

Temporary connection during placement

Radiopaque marker

Supports visualisation under imaging

Compatible Urotech ureteral stents

○ WHITE-STAR® ● YELLOW-STAR® ● GREEN-STAR® ● TUMORSTENT®



Ordering information

UROCKIT® Steering Device

Order code	Compatibility	Working length (cm)	Wire diameter (inch)
SD-0750-M	6, 7 und 8 CH/FR Ureterschienen	42	0.035
SD-0550-M	5 CH/FR Ureterschienen	42	0.035



UROCKIT® compatible ureteral stents

The focus products shown represent a selection of compatible ureteral stents for use with item SD-0750-M. Further combinations of the UROCKIT® Steering Device and compatible ureteral stents are available on request. Please contact your responsible field sales representative.

 WHITE-STAR®	 YELLOW-STAR®	 GREEN-STAR®	 TUMORSTENT®
Variants without wire	Variants without wire	Variants without wire	Variants without wire
EN-410624-F	EP-410624	EG-410624	TU-410624
EN-410626-F	EP-410626	EG-410626	TU-410626
EN-410628-F	EP-410628	EG-410628	TU-410628
EN-410630-F	EP-410630	EG-410630	TU-410724
EN-410724-F	EP-410724	EG-410724	TU-410726
EN-410726-F	EP-410726	EG-410726	TU-410728
EN-410728-F	EP-410728	EG-410728	TU-410730
EN-410730-F	EP-410730	EG-410730	TU-410824
EN-410824-F	EP-410824	EN-410826	TU-410826
EN-410826-F	EP-410826	EN-410828	TU-410828
EN-410828-F	EP-410828		TU-410830
EN-410830-F	EP-410830		
Variants with wire	Variants with wire		
EN-360624-F	EP-360624		
EN-360626-F	EP-360626		
EN-360628-F	EP-360628		
EN-360630-F	EP-360630		
EN-360724-F	EP-360724		
EN-360726-F	EP-360726		
EN-360728-F	EP-360728		
EN-360730-F	EP-360730		
EN-360824-F	EP-360824		
EN-360826-F	EP-360826		
EN-360828-F	EP-360828		
EN-360830-F	EP-360830		

Article number key

EN-410624-F

Product family ——— Variant with thread
Diameter (CH/FR) ——— Shaft length (cm)



Urotech GmbH
Medi-Globe-Str. 1-5
83101 Rohrdorf, Germany
Tel: +49 (0) 8032 - 973-333
Fax: +49 (0) 8032 - 973-211
E-Mail: sales@urotech.com
www.urotech.de



eIFU